



ZHANGJIAGANG ANNECY MACHINERY CO., LTD.

CATALOGUE & QUOTATION

Images	Description	Model No	Details
	Baby bed	AC-BB001	(1) Stainless steel 201 bed frame. (2) With 4 silent castor. (3) With mattress. (4) Size:750*450*950mm (5) N.W.: 11kgs (6) G.W.: 13kgs (7) MOQ: 3 pieces (8) Production ability: 4500 pieces/month (9) With 4 castors, and mattress.
	Baby Bed	AC-BB002	Specification: L820*W460*H750/1000mm Material: full stainless steel with transparent GRP and castor Optional accessory: Mattress
	Baby Bed	AC-BB003	Specification: L820*W460*H650cm
	Baby Bed	AC-BB004	Specification: L820*W460*H650/900mm full stainless steel with transparent GRP and castor Optional accessory: Mattress
	Baby Bed	AC-BB005	Specification: L730*W410*H610mm full stainless steel with transparent GRP and castor Optional accessory: Mattress

	Baby Bed	AC-BB006	<p>Features:</p> <p>Products can be tilted left and right, in the range of inclination can be any angle positioning, tilt range: $\pm 15^\circ$</p> <p>Low—High:780—980 mm , the height can be adjusted ,easy for mother and baby to communicate.</p> <p>Size: L905*W535*H845mm</p> <p>Groove in the Chassis guard temporary display the baby wash</p> <p>Baby basin can be put in anyway for cleaning and refresh.</p> <p>The castor with brake, without brake assembly opposite angle, it's much more stable and reliable .</p> <p>This product is mainly used to bathe the baby, can also be used for leisure, sleep.</p>
	Baby Bed	AC-BB007	<p>Height adjustable by hydraulic system , 1045mm-855mm</p> <p>.Tilt angle 15° ,easy for mother and baby to communicate and feed</p> <p>Humanity functions :</p> <ol style="list-style-type: none"> 1. Music playing : With a SD memory card ,could storage music in MP3 format. 2. Weight scale : Measure weight at anytime 3. Atmosphere temperature measuring .
	Baby Bed	AC-BB008	<p>Specification: L1850*W750*H550mm</p> <ol style="list-style-type: none"> 1.Rectangular frame of Head/Foot board with crank 2. Stamping bed surface 3. Up-down lift side rail
	Baby Bed	AC-BB009	<p>Outside size L1850*W750*H550mm</p> <p>Inside size of mattress base L1850*W650mm</p> <p>bedboard Steel with punching holes</p> <p>Head & foot board Rectangular frame</p> <p>Guardrails Up-down lift side rail</p> <p>Crank one piece</p> <p>Wheels 5" castor with brakes</p> <p>Safe working load The quality attains national standard.It is equipped with wheel which is quiet and noiseless and able to support over 160kgs.</p>
	Baby Bed	AC-BB010	<p>Aluminum alloy head & foot boards with laminated panels, can be removable</p> <ul style="list-style-type: none"> • Aluminum alloy side rail, 5-staff, collapsed type with safe lock • Bed sheet top made from perforated steel,epoxy powder coated and baked finish • Bed framework made from mild steel ERW rectangular tubes • Dimension: 1800*780*500mm • 2 functions: <p>Backrest lifting angle: $80^\circ \pm 5^\circ$</p> <p>Kneerest lifting angle: $45^\circ \pm 5^\circ$</p> <ul style="list-style-type: none"> • With shoes holder • Gross weight:105kg.
	Baby Bed	AC-BB011	<ol style="list-style-type: none"> (1) Stainless steel 201bed frame. (2) With 4 silent castor. (3) With mattress. (4) Size:750*450*950mm (5) N.W.: 11kgs (6) G.W.: 13kgs (7) MOQ: 3 pieces (8) Production ability: 4500 pieces/month (9) With 4 castors, and mattress.

	Baby Bed	AC-BB012	Product Size: 750*450*950mm Made by stainless steel 201 Stainless steel baby bed Full stainless steel bassinet and baby cot With 4 castors, and mattress.
	Baby Bed	AC-BB013	Size: 800*500*(780-980)mm Painted steel with basket Small castors with brakes
	Baby Bed	AC-BB014	Size: 800*500*780mm Powder coated steel structure Small castors with brakes
	Baby Bed	AC-BB015	Wooden material baby cot Easy assemble, with one big drawer Four small castors noiseless and with brakes Size: 810*520*710mm
	Baby Bed	AC-BB016	Stainless steel baby cot Size: 730*460*1030mm Mattress: high density sponge + PU cover Colour: blue, green, black Weight: 25 kgs
	Baby Bed	AC-BB017	Size: 1210*690*600mm Flat bed Structure: Wholly molded cold steel sleeping platform Electrophoresis and power coating Head & Foot board : Steel Side rails: full length protection, steel, sliding up and down Castors: 125mm castors, two with brakes Loading capacity: 160kgs