






**ZHANGJIAGANG ANNECY MACHINERY CO., LTD.**

**Emergency trolley Catalogue & Quotation**


Images	Model No	Details	Size
	AC-ET001	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole	850*520*960mm
	AC-ET002	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.C.P.R board 9.Oxygen Tank Holder	750*475*930mm






	AC-ET003	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers and two doors 2.3*3 segregation board in every drawer 3.I.V. pole	850*520*960mm
	AC-ET004	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.C.P.R board 9.Oxygen Tank Holder	750*475*930mm




	AC-ET005	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 626 * 475*930mm Packing size : 675*520*970mm Gross weight: 36kg
	AC-ET006	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder	Specification: 650 * 480*930mm Packing size : 720*540*970mm Gross weight: 41kg




	AC-ET007	Mainly made of plastic, aluminum and steel, ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 750 * 475*930mm Packing size : 770*540*940mm Gross weight: 37.75KG
	AC-ET008	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder	Specification: 850 * 520*930mm Packing size : 890*570*980mm Gross weight: 42KG

	AC-ET009	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 750 * 475*930mm Packing size : 730*530*940mm Gross weight: 37.75
	AC-ET010	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf	Specification: 850* 520*930mm Packing size : 770*540*940mm Gross weight: 42
	AC-ET011	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 750 * 475*930mm Packing size : 730*530*940mm Gross weight: 37.75
	AC-ET012	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.C.P.R board 9.Oxygen Tank Holder	Specification: 850 * 520*930mm Packing size : 770*540*940mm Gross weight: 42
	AC-ET013	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two small drawers one middle size drawer and one tray 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 750 * 475*930mm Packing size : 730*530*940mm Gross weight: 37.75
	AC-ET014	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder	Specification: 850 * 520*930mm Packing size : 770*540*940mm Gross weight: 42














	AC-ET015	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 625 * 475*930mm Packing size : 675*520*970mm Gross weight: 36
	AC-ET016	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder	Specification: 650 * 480*930mm Packing size : 720*540*970mm Gross weight: 41
	AC-ET017	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers . 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 626* 475*930mm Packing size : 675*520*970mm Gross weight: 36
	AC-ET018	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder	Specification: 650 * 480*930mm Packing size : 720*540*970mm Gross weight: 41
	AC-ET019	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers. 2.3*3 segregation board in every drawer 3.I.V. pole	Specification: 626 * 475*930mm Packing size : 675*520*970mm Gross weight: 36kg
	AC-ET020	4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder	Specification: 650 * 480*930mm Packing size : 720*540*970mm Gross weight: 41kg









	AC-ET021	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3'3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder	850*520*950mm
	AC-ET022		750*475*930mm



	AC-ET023	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3'3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.Defibrillator shelf	850*520*950mm
	AC-ET024		750*475*930mm
	AC-ET025	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3'3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder 11.Needle Disposal Holder	625*475*920mm

	<p>AC-ET026</p>	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	<p>625*475*920mm</p>
	<p>AC-ET027</p>	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	<p>625*475*920mm</p>
	<p>AC-ET028</p>	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder 11.Needle Disposal Holder</p>	<p>673*480*930mm</p>
	<p>AC-ET029</p>	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	<p>630*470*930mm</p>
	<p>AC-ET030</p>	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	<p>673*480*930mm</p>
	<p>AC-ET031</p>	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Double trash cans 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	<p>630*470*930mm</p>
			
	<p>AC-ET032</p>	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Single waste pail 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	<p>750*475*920mm</p>
			

	AC-ET033	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Single waste pail 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	750*475*920mm
	AC-ET034	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Single waste pail 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	750*470*930mm
			
	AC-ET035	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Two drawers one door 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Storage box 7.Single waste pail 8.Defibrillator shelf 9.C.P.R board 10.Oxygen Tank Holder</p>	750*470*930mm
	AC-ET036	<p>Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Five drawers 2.3*3 segregation board in every drawer 3.I.V. pole 4.Power supply equipped 5.Sliding Side Shelf working table 6.Side tumbling type material box 7.Double trash cans 8.C.P.R board 9.Oxygen Tank Holder</p>	<p>Specification: 750 * 475*930mm Packing size : 730*530*940mm Gross weight: 37.75</p>
	AC-ET037		<p>Specification: 850 * 520*930mm Packing size : 770*540*940mm Gross weight: 42</p>