












ZHANGJIAGANG ANNECY MACHINERY CO., LTD.

CATALOGUE & QUOTATION

Images	Description	Model No	Details
	Foldable stretcher	AC-FS001	1) Product size (L x W x H): 185 x 50x 20cm 2) Folded size (L x W x H): 92.5x 50 x 9cm 3) Packing size (2pcs/carton): 96 x 54 x 16cm 4) N.W.: 6kg G.W.:7 kg 5) Load bearing: less than 159kg
	Foldable stretcher	AC-FS002	1) Product size (L x W x H): 185 x 50 x 18cm 2) Folded size (L x W x H): 92.5x 50 x 9cm 3) Packing size (2pcs/carton): 96 x 54 x 16cm 4) N.W.: 6kg G.W.: 7kg 5) Load bearing: less than 159kg
	Foldable stretcher	AC-FS003	1) Product size (L x W x H): 185 x 50 x 24cm 2) Folding size (L x W x H): 92.5 x 50 x 11cm 3) Packing size (1pcs): 96 x 54 x 12cm
	Foldable stretcher	AC-FS004	1) Product size (L x W x H):185 x 50x 3cm 2) Folded size (L x W x H):92.5 x 50 x 7cm 3) Packing size(2pcs/carton): 96 x 54 x 16cm 4) N.W.: 4kg 5) G.W.: 5.1kg 6) Load bearing: less than 159kg
	Foldable stretcher	AC-FS005	1) Product size (L x W x H):185 x 50 x 3cm 2) Folding size (L x W x H):92.5 x 50 x 7cm 3) Packing size (2pcs/carton): 96 x 54 x 16cm 4) N.W.: 6.5kg G.W.: 7.6kg 5) Load bearing: less than 159kg
	Foldable stretcher	AC-FS006	1) Product size (L x W x H):185 x 50 x 3cm 2) Folded size (L x W x H):92.5 x 50 x 10cm 3) Packing size(2pcs/carton): 96 x 54 x 16cm 4) N.W.: 6.5kg 5) G.W.: 7.6kg 6) Load bearing: less than 159kg, steel
	Foldable stretcher	AC-FS007	1) Product size (L x W x H): 208 x 54.5 x 13cm 2) Folding size (L x W x H): 104 x 9 x 17cm 3) Packing size (2pcs): 108 x 18 x 18.5cm 4) N.W.: 4.7kg 5) G.W.: 5.3kg 6) Load bearing: less than 159kg

	Foldable stretcher	AC-FS008	<ol style="list-style-type: none"> 1) Product size (L x W x H): 214 x 54.5 x 12cm 2) Folded size (L x W x H): 58 x 21 x 14cm 3) Packing size(1pc/carton): 67 x 30 x17cm 4) N.W.: 5kg 5) G.W.: 6.6kg 6) Load bearing: less than 159kg
	Foldable stretcher	AC-FS009	<ol style="list-style-type: none"> 1) Product size (L x W x H): 215x 54.5 x 13cm 2) Folded size (L x W x H): 182x 17 x 9cm 3) Packing size (1pc/carton):185x 19 x 15cm 4) N.W.: 6.5kg 5) G.W.: 7.8kg 6) Load bearing: less than 159kg
	Foldable stretcher	AC-FS010	<ol style="list-style-type: none"> 1) Product size (L x W x H): 208x 54.5 x 15cm 2) Folded size (L x W x H): 104x 17 x 9cm 3) Packing size (1pc/carton): 18.5x 18 x 108cm 4) N.W.: 5kg 5) G.W.: 5.6kg 6) Load bearing: less than 159kg
	Foldable stretcher	AC-FS011	<ol style="list-style-type: none"> 1) Product size (L x W x H): 215x 54.5 x 15cm 2) Folded size (L x W x H): 26x 17 x 55cm 3) Packing size (1pc/carton):30x 17 x 58cm 4) N.W.: 7.7kg
	Foldable stretcher	AC-FS012	<p>Product Size: 210*52*14cm Packing Size: 100*52*14cm Seat-H: 52cm Folded into wheelchair Back-H: 137cm Load Bearing: 159Kg N.W.: 8KG G.W.: 9.5KG</p>
	Foldable stretcher	AC-FS013	<ol style="list-style-type: none"> 1) Product size (L x W x H): 220 x 55 x 15cm 2) Folded size (L x W x H): 180 x 12 x 15cm 3) Packing size (1pc): 183 x 14 x 16cm 4) N.W.: 7.5kg G.W.: 9.5kg 5) Load bearing: less than 159kg
	Foldable stretcher	AC-FS014	<ol style="list-style-type: none"> 1) Product size (L x W x H): 220 x 55 x 15cm 2) Folded size (L x W x H): 95 x 20 x 9cm 3) Packing size (1pc/bag): 98 x 22x 11cm 4) N.W.: 7.5kg G.W.: 8.2kg 5) Load bearing: less than 159kg
	Foldable stretcher	AC-FS015	<ol style="list-style-type: none"> 1) Product size (L x W x H): 185 x 48 x 3cm 2) Folded size (L x W x H): 46 x 48 x 14cm 3) Packing size (2pcs/carton): 50x 48 x 15cm 4) N.W.: 4.5kg G.W.: 6.2kg 5) Load bearing: less than 159kg

	Foldable stretcher	AC-FS016	<p>1) Product size (L x W x H): 185 x 48 x 20cm 2) Folding size (L x W x H): 46 x 48 x 14cm 3) Packing size (2pcs/carton): 50x 48x 15cm 4) N.W.: 7kg G.W.: 8.2kg 5) Load bearing: less than 159kg 6) G.W.: 8.4kg 6) Load bearing: less than 159kg</p>
	Foldable stretcher	AC-FS017	<p>1)Product size (L x W x H):208x53x35cm 2)Folded size:108x53x12cm 3)N.W.: 8.4kg G.W.: 9KG 4)Packing size (2pcs/carton): 112x55x22cm 5)Load bearing: less than 159kg</p>
	Foldable stretcher	AC-FS018	<p>This foldable stretcher with wheels, outriggers and can be folded. Foldable stretcher is made of high strength aluminum alloy material and the waterproof leather materials. Foldable stretcher has characterizes of light-weighted, small-sized, easy-carried, reliable quality, use-safely and easy for sterilization, etc. Foldable stretcher is commonly used in funeral house to carry and transport the dead bodies.</p> <p>Technical Features: 1) Product size (L x W x H): 185 x 50x 20cm 2) Folded size (L x W x H): 92.5x 50 x 12cm 3) Packing size (1pc/carton): 96 x54x 20cm 4) N.W.: 8kg G.W.:9.5 kg 5) Load bearing: less than 159kg</p>
	Foldable stretcher	AC-FS019	<p>1) Product size (L x W x H): 220 x 55 x Open size: 230*17*55cm Folding size: 103*14*22cm Packing size: 105*16*24cm Loading Capacity: 159KG N.W.: 9.4KG G.W.: 10.5KG</p> <p>Stretcher used for seriously injured person ,it is made of aluminum alloy ,it can save the rescue time ,it has characteristics of durable water -proof ,small-size as it can piled up ,so it can save he storage space, what more is that it can avoid the pollution of the flood</p>
	Foldable stretcher	AC-FS020	<p>1)Product size (L x W x H):187*50*18 cm 2)packing size:190*52*20cm 3)N.W.: 5kg G.W.: 5.5kg 4)Load bearing: less than 159kg</p> <p>Stretcher used for seriously injured person ,it is made of aluminum alloy ,it can save the rescue time ,it has characteristics of durable water -proof ,small-size as it can piled up ,so it can save he storage space, what more is that it can avoid the pollution of the flood</p>
	Foldable stretcher	AC-FS021	<p>1)Product size (L x W x H):202x53x9.5cm 2)Folded size:100x53x13cm 3)N.W.: 10kg G.W.: 10.5kg 4)Packing size (2pcs/carton): 112x55x22cm 5)Load bearing: less than 159kg</p>
	Foldable stretcher	AC-FS022	<p>Product Size: 230*55*13.5cm Folded size: 100*14*12cm Packing Size: 105*16*14cm Loading bearing: <=159kgs N.W.:9.5KG G.W.: 10KG</p> <p>Detailed Product Description: *It is made of high strength aluminum alloy material and made by PU leather materials. *This stretcher can be folded by lengthwise with flexible hand.r. *It has characterizes of light-weighted, small-sized, easy-carried, reliable quality, use-safely and easy for sterilization, etc. *It is commonly used in hospital, gymnasium, ambulance, and battle fields to carry patients and wounded person.</p>

	Soft Stretcher	AC-FS023	<p>Soft stretcher can be making according to your requirement. Unfolding size:187*42cm. Self weight:7kg. Bearing weight:less than 204kg Net weight:8kg Gross weight:9kg. Usually it is not with bag,if you need ,we can customized.</p>
	Soft Stretcher	AC-FS024	<p>Soft stretcher can be making according to your requirement. Unfolding size:212*75cm. Folding size:88*24cm Self weight:7kg. Bearing weight:less than 159kg Net weight:9kg Gross weight:10kg. It is with a bag</p>
	Vacuum Stretcher	AC-FS025	<p>Product Description:</p> <p>TPU material: Tough: soft as silk skin friendly material on the patient side. Extremely puncture resistant on the outside.</p> <p>Valve automatically closes when suction stops, no risk of losing vacuum if slipping with the vacuum pump. It can be operated with gloves. To evacuate: turn valve to the right. To release: turn valve to the left.</p> <p>The patented Multi Chamber System which keeps to beads in the right position and makes fixation quick, safe and easy.</p> <p>Unique handle and strapping system: color coded straps which are adjustable from both sides; Heavy duty cushioned handles provides optimal weight distribution and are easy on the hands; Straps and handles can easily be removed for repositioning, cleaning or replacement.</p> <p>It is X-RAY translucent and can be used for MRI. Body Shaped Vacuum Mattress (Adult) Width: 88cm Length: 200cm Weight: ~7.3kg</p>