











ZHANGJIAGANG ANNECY MACHINERY CO.LTD

CATALOGUE & QUOTATION

Images	Description	Model No
	Attendant chair	AC-AC001 Specification: Frame is made steel tube with plastic coating. Mattress is made of PVC artificial leather. The Attendant chair can be open as a simple hospital bed, it can be folded into chair. This chair is for the person who accompany patients in hospital, it can be fold as a chair in day time, while unfold as a bed in night.Size Folding:620*720*940mm,unfolding: 620*1900*420mm Packing size 650*770*1380mm./2 pcs Material Steel powder coated and artificial leather CBM 0.7/2PC, Gross weight 26kg/pc, 520kg/20pc Accessories One pillow; Three sections; Handrails with ABS Engineering plastic arm rest; Leg section with four castors
	Attendant chair	AC-AC002 Unfolded size:2060*650*450mm,Folding size:L650*W780*H1020mm,material:Aluminium alloy armrest,mute wheels,Powder costed steel frame waterproof PVC cover with foam inside.wheels:6pcs.NW:38KG,Packing size:130*66*38cm.
	Attendant chair	AC-AC003 Specification: 600*720/1900*250/380mm
	Transfusion chair	AC-TC001 Specification : L680*W800*H1100mm
	Transfusion chair	AC-TC002 Specification: L730*W780*H1050mm

	Waiting chair	AC-WTC001	<p>Size 181*70*79cm</p> <p>Specification Seat board is made of high-quality cold rolling steel plate</p> <p>Legs, arms, and the edges of backrest are all chromed</p> <p>With three seats</p> <p>gw:32kg,volume:0.1cbm. 3 pcs in one carton</p>
	Waiting chair	AC-WTC002	<p>Size 175*66*78cm</p> <p>Specification Seat board is made of high-quality cold rolling steel plate</p> <p>Legs, arms, and the edges of backrest are all chromed and with PU cover</p> <p>With three seats two pcs in one carton,GW:36KG.volume:0.09cbm</p>
	Nursing stool	AC-NS001	<p>Specification : H480~650mm by Gas Spring,Stool suface:φ 300mm.</p> <p>Description: the fabric is pressed PVC artificial leather and steel frame chromed</p>
	Nursing stool	AC-NS002	<p>Specification : H480~650mm by Gas Spring,Stool suface:φ 300mm.</p> <p>Description: the fabric is pressed PVC artificial leather and steel frame chromed</p>
	Nursing stool	AC-NS003	<p>Specification : H480~650mm by Gas Spring,Stool suface:φ 300mm.</p> <p>Description: the fabric is pressed PVC artificial leather and steel frame chromed</p>



	Nursing stool	AC-NS004	<p>Seat:pu leather + high resilient foam+abs board ,size:435*400mm Back:material:pu leather + high resilient foam+abs board ,size:400*240mm.</p> <p>Height: 470-590mm. 485-605mm.</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel,aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS005	<p>Seat:pu leather + high resilient foam+abs board ,size:435*400mm Back:material:microfiber leather + high resilient foam+abs board ,size:380*310mm.</p> <p>Height: 470-590mm. 485-605mm.</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel,aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS006	<p>Seat:pu leather + high resilient foam+abs board ,size:450*410mm Back:material:pu leather + high resilient foam+abs board ,size:400*240mm.</p> <p>Height: 450-570mm. 465-585mm.</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel,aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS007	<p>Seat:Microfiber leather + high resilient foam+abs board ,size:460*420mm</p> <p>Back:material:Microfiber leather + high resilient foam+abs board ,size:410*250mm.</p> <p>Height: 460-580mm. 475-595mm.</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel,aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS008	<p>Seat:pu leather + high resilient foam+abs board ,size:450*410mm</p> <p>Back:material:pu leather + high resilient foam+abs board ,size:400*240mm.</p> <p>Height: 460-580mm</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, radius:280mm</p> <p>Castor:pu mute. 50mm</p>

	Nursing stool	AC-NS009	Seat:pu leather + high resilient foam+abs board ,size:450*410mm Back:material:pu leather + high resilient foam+abs board ,size:370*300mm. Height: 460-580mm Travel distance:120mm, Surface treatment:chromed Material:chromed steel, radius:280mm Castor:pu mute. 50mm
	Nursing stool	AC-NS010	Seat:Microfiber leather + high resilient foam+abs board ,size:460*420mm Back:material:Microfiber leather + high resilient foam+abs board ,size:400*250mm. Height: 485-605mm Travel distance:120mm, Surface treatment:chromed Material:chromed steel, radius:280mm Castor:pu mute. 50mm
	Nursing stool	AC-NS011	Seat:Microfiber leather + high resilient foam+abs board ,size:460*420mm Back:material:Microfiber leather + high resilient foam+abs board ,size:380*310mm. Height: 485-605mm Travel distance:120mm, Surface treatment:chromed Material:chromed steel, radius:280mm Castor:pu mute. 50mm
	Nursing stool	AC-NS012	Seat:pu leather + high resilient foam ,size:440*400mm Back:material:pu leather + high resilient foam ,size:370*270mm. Height: 420-510mm, 435-555mm Travel distance:120mm, Surface treatment:chromed Material:chromed steel, radius:280mm Castor:pu mute. 50mm
	Nursing stool	AC-NS013	Seat:Microfiber leather + high resilient foam ,size:450*410mm Back:material:Microfiber leather + high resilient foam ,size:380*280mm. Height: 430-550mm, 445-565mm Travel distance:120mm, Surface treatment:chromed Material:chromed steel, radius:280mm Castor:pu mute. 50mm

	Nursing stool	AC-NS014	<p>Seat:pu leather + high resilient foam+abs board, size:450*410mm</p> <p>Back:material:pu leather + high resilient foam+abs board, size:370*300mm.</p> <p>Height: 450-570mm, 465-585mm</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS015	<p>Seat:Microfiber leather + high resilient foam+abs board , size:460*420mm</p> <p>Back:material:Microfiber leather + high resilient foam+abs board , size:380*310mm.</p> <p>Height: 460-580mm. 475-595mm.</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS016	<p>Seat:pu leather + high resilient foam+abs board , size:Φ 360mm</p> <p>Back:material:self-skin pu + high resilient foam+abs board , size:450*50mm.</p> <p>Height: 460-570mm.</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS017	<p>Seat:pu leather + high resilient foam+abs board , size:Φ 360mm</p> <p>Back:material:self-skin pu + high resilient foam+abs board , size:450*80mm.</p> <p>Height: 450-570mm. 480-600mm</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS018	<p>Seat:pu leather + high resilient foam+abs board , size:Φ 360mm</p> <p>Back:material:self-skin pu + high resilient foam+abs board , size:290*110mm.</p> <p>Height: 440*560mm, 470-590mm</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>

	Nursing stool	AC-NS019	<p>Seat:pu leather + high resilient foam+abs board ,size:Φ 360mm</p> <p>Back:material:self-skin pu + high resilient foam+abs board ,size:450*50mm.</p> <p>Height:570-770mm.</p> <p>Footring:aluminium, chromed steel, Φ 400mm</p> <p>Travel distance:200mm, Surface treatment:chromed</p> <p>Material:chromed steel, aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS020	<p>Seat:pu leather + high resilient foam+abs board ,size:Φ 330mm</p> <p>Height:440-560mm.</p> <p>Footring:aluminium, chromed steel, Φ 400mm</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, radius:220mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS021	<p>Seat:pu leather + high resilient foam+abs board ,size:Φ 330mm</p> <p>Height:420-540mm.</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, radius:240mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS022	<p>Seat:self-skin pu , size:340*230mm</p> <p>Height:440-580mm. back size:150*150mm</p> <p>Footring:aluminium, chromed steel, Φ 400mm</p> <p>Travel distance:120mm, Surface treatment:chromed</p> <p>Material:chromed steel, aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS023	<p>Seat:self-skin pu , size:415*390mm</p> <p>Height:450-590mm. back size:390*255mm</p> <p>Footring:aluminium, chromed steel, Φ 400mm</p> <p>Travel distance:140mm, Surface treatment:chromed</p> <p>Material:chromed steel, aluminium, radius:280mm</p> <p>Castor:pu mute. 50mm</p>

	Nursing stool	AC-NS024	<p>Seat:self-skin pu , size:415*390mm Height:450-590mm. back size:370*220mm Footring:aluminium, chromed steel, Φ400mm Travel distance:140mm, Surface treatment:chromed Material:chromed steel, aluminium, radius:280mm Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS025	<p>Seat:self-skin pu+abs board , size:420*350mm Height:460-600mm. Footring:aluminium, chromed steel, Φ400mm Travel distance:140mm, Surface treatment:chromed Material:chromed steel, aluminium, radius:280mm Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS026	<p>Seat:self-skin pu+abs board , size:420*350mm Height:460-600mm. Footring:aluminium, chromed steel, Φ400mm Travel distance:140mm, Surface treatment:chromed Material:chromed steel, aluminium, radius:280mm Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS027	<p>Seat:self-skin pu+abs board , size: Φ320mm Height:425-565mm. Footring:aluminium, chromed steel, Φ400mm Travel distance:140mm, Surface treatment:chromed Material:chromed steel, aluminium, radius:280mm Castor:pu mute. 50mm</p>
	Nursing stool	AC-NS028	<p>Seat:self-skin pu , size:460*290mm Height:420-560mm. back size:380*220mm Footring:aluminium, chromed steel, Φ400mm Travel distance:140mm, Surface treatment:chromed Material:chromed steel, aluminium, radius:280mm Castor:pu mute. 50mm</p>

	Nursing stool	AC-NS029	Seat:self-skin pu , size:480*450mm Height:450-590mm. back size:440*300mm Footring:aluminium, chromed steel, Φ400mm Travel distance:140mm, Surface treatment:chromed Material:chromed steel, aluminium, radius:280mm Castor:pu mute. 50mm
	Nursing stool	AC-NS030	Seat:self-skin pu , size:350*310mm Height:45-575mm. back size:285*160mm Footring:aluminium, chromed steel, Φ400mm Travel distance:140mm, Surface treatment:chromed Material:chromed steel, aluminium, radius:280mm Castor:pu mute. 50mm