
















ZHANGJIAGANG ANNECY MACHINERY CO., LTD.

LUXURIOUS TROLLEY CATALOGUE & QUOTATION

Images	Description	Model No	Details
	Emergency Trolley	AC-LET01	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer 8.Calm and good rolling 4" casters, two casters with brake 9.Size:615*550*1010mm <p>Standard accessories: Defibrillator shelf 1 pcs, IV pole, 1 hanging box, 2 dust box without cover, ABS CPR board Oxygen tank holder, 220V power outlet with hooks.</p> 
	Emergency Trolley	AC-LET03	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer 8.Calm and good rolling 4" casters, two casters with brake 9.Two needle boxes 10.Size:615*550*1010mm <p>Standard accessories: Defibrillator shelf 1 pcs, IV pole, 1 hanging box, 2 dust box without cover, ABS CPR board Oxygen tank holder, 220V power outlet with hooks.</p>
	Emergency Trolley	AC-LET04	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer 8.Calm and good rolling 4" casters, two casters with brake 9.Size:615*550*1010mm <p>Standard accessories: Defibrillator shelf 1 pcs, IV pole, 1 hanging box, 2 dust box without cover, ABS CPR board Oxygen tank holder, 220V power outlet with hooks.</p>

	<p>Medicine Trolley</p>	<p>AC-LMT01</p>	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer 8.Calm and good rolling 4" casters, two casters with brake 9.Two needle boxes 10.Size:615*550*1010mm <p>Standard accessories: 1 hanging box, 2 dust box without cover</p>
			
	<p>Medicine Trolley</p>	<p>AC-LMT02</p>	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer and one side door 8.Calm and good rolling 4" casters, two casters with brake 9.Two needle boxes 10.Size:920*640*1060mm <p>Standard accessories: 3 hanging box, 2 dust box without cover</p>
	<p>Medicine Trolley</p>	<p>AC-LMT03</p>	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer 8.Calm and good rolling 4" casters, two casters with brake 9.Two needle boxes 10.Size:550*470*1010mm <p>Standard accessories: 1 hanging box, 2 dust box without cover</p>

	Anesthesia Trolley	AC-LAT01	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer 8.Calm and good rolling 4" casters, two casters with brake 9.Two needle boxes 10.Eight Multi Bin Container 11.Size:615*550*1470mm <p>Standard accessories:</p> <p>1 hanging box, 2 dust box without cover</p>
	Anesthesia Trolley	AC-LAT03	<p>Technical Parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating , all side made of ABS molding material 3.Anti-collision corner each side 4.3 ABS utility containers on left side storage 5.Right side 1 hanging box and 2 dust box 6.Centralized locking 7.2 small drawers, 2 middle drawer and 1 large drawer 8.Calm and good rolling 4" casters, two casters with brake 9.Two needle boxes 10.Ten Multi Bin Container 11.Size:615*550*1470mm <p>Standard accessories:</p> <p>1 hanging box, 2 dust box without cover</p>
	Emergency Trolley	AC-LET02	<p>Technical parameter:</p> <ol style="list-style-type: none"> 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating, All side made of ABS molding material 3.Anti collision corner each side 4.Centralized locking 5.2 small drawers, 2 middle drawer and 1 large drawer 6.Calm and good rolling 4" castors, two castors with brake 7.Size: 810*520*960mm <p>Standard accessories:</p> <p>Defibrillator shelf 1 pcs, IV pole, 1 hanging box , 2 dust box without cover, ABS CPR board Oxygen tank holder, 220V power outlet with hooks.</p>

	<p>Medicine Trolley</p>	<p>AC-LMT04</p>	<p>Technical parameter: 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating, All side made of ABS molding material 3.Anti collision corner each side 4.Centralized locking 5.2 small drawers, 2 middle drawer and 1 large drawer 6.Calm and good rolling 4" castors, two castors with brake 7.Size: 810*520*960mm Standard accessories: 1 hanging box , 2 dust box without cover</p>
	<p>Anesthesia Trolley</p>	<p>AC-LAT02</p>	<p>Technical parameter: 1.ABS molding working top frame 2.1.2 Cold frame with Electrostatic powder coating, All side made of ABS molding material 3.Anti collision corner each side 4.Centralized locking 5.2 small drawers, 2 middle drawer and 1 large drawer 6.Calm and good rolling 4" castors, two castors with brake 7.Eight Multi Bin Container 8.Size: 810*520*1420mm Standard accessories: 1 hanging box , 2 dust box without cover</p>
			