






















ZHANGJIAGANG ANNECY MACHINERY CO., LTD.



CATALOGUE & QUOTATION


Images	Description	Model	Details
	five functions hospital beds for sale	AC-MB001	<p>Specification: L2150*W950*H400-720mm</p> <ol style="list-style-type: none"> 1. Detachable head/foot board 2. Collapsible side rails by steel coated 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	five functions hospital beds for sale	AC-MB002	<p>Specification: L2150*W950*H400-720mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	five functions hospital beds for sale	AC-MB003	<p>Specification:L2150*W900*H450/750mm</p> <ol style="list-style-type: none"> 1.PP material collapsible side rails 2.Detachable ABS head/foot board 3.I.V. Rod & Drainage hook 4.5"Castor with brakes 5.Weight capacity:over 250kgs <p>Function:</p> <ol style="list-style-type: none"> 1) Backrest max upward angle: 75° 2) Footrest max upward angle: 45° 3) Height adjustment: 450-700mm 4)Trendelenburg:0-10° 5)Reversed Trendelenburg:0-10°
	three functions medical bed for sale	AC-MB004	<p>Specification:L2120*W900*H450/750mm</p> <p>Back Rest 0-75°</p> <p>Leg Rest 0-45°</p> <p>Maximum safe load 250KG</p> <p>Diameter of Centrally Control Silent and Dustproof Castor 100mm</p> <p>Special Allocation</p> <ul style="list-style-type: none"> • I.V.Pole (optional) • Drainage Hook • Angle Scale • Mattress (optional) • Detachable PP Head and foot board • Aluminum Alloy Side Rail • Top Coat: Electrostatic Powder coating • Bed frame: Cold-Rolled Steel Tube • Bed base: Cold-Rolled Steel Sheet
	three functions medical bed for sale	AC-MB005	<p>Specification: L2150*W1050*H400-720mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. Center control lock 6. Optional accessory: Mattress & over-bed table & bedside cabinet





	<p>three functions medical bed for sale</p>	<p>AC-MB006</p>	<p>Specification: L2150*W1050*H460-720mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible ABS guardrails 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	<p>three functions medical bed for sale</p>	<p>AC-MB007</p>	<p>Specification: L2150*W1050*H460-720mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible aluminum side rails with diner table 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. Center control lock 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	<p>three functions hospital bed price</p>	<p>AC-MB008</p>	<p>Specification: L2150*W1050*H460-720mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by ABS with diner table 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. Center control lock 6. Optional accessory: Mattress & over-bed table & bedside cabinet

	three functions hospital bed price	AC-MB009	<p>Specification: L2150*W950*H350-690mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	three functions hospital bed price	AC-MB010	<p>Specification: L2150*W950*H350-690mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	two functions patient bed	AC-MB011	<p>Specification: L2150*W950*H500mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum with diner table 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with two brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	two functions patient bed	AC-MB012	<p>Specification: L2150*W950*H500mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum with diner table 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with two brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet
	two functions patient bed	AC-MB013	<p>Specification: L2150*W950*H500mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. Center control lock 6. Optional accessory: Mattress & over-bed table & bedside cabinet <p>ABS double crank hospital bed abs double crank hospital bed abs double crank hospital bed Mattress, Over-bed table, Beside cabinet</p>

	two functions patient bed	AC-MB014	<p>Specification: L2150*W950*H500mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with two brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet <p>Back-rest lifting 0-75°(±2%) Knee-rest lifting 0-45°(±2%) Optional Accessory: ABS double crank hospital bed abs double crank hospital bed abs double crank hospital bed Mattress, Over-bed table, Beside cabinet</p>
	two functions hospital patient bed	AC-MB015	<p>Specification: L2150*W950*H500mm</p> <p>Standard Accessory: Detachable ABS head/foot board Collapsible aluminum alloy guardrails Diner table</p> <p>Integrated stamping frame and mattress base are made of steel epoxy I.V Pole Drainage hook with PU cover Heavy duty caster with brakes 5 inch Manual crank system 2set Optional accessory: Mattress over-bed table bedside cabinet</p>
	two functions hospital patient bed	AC-MB016	<p>Specification: L2150*W950*H500mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 2. Collapsible side rails by aluminum 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with two brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet <p>Back-rest lifting 0-75°(±2%) Knee-rest lifting 0-45°(±2%) Optional Accessory: second hand hospital bed second hand hospital bed second hand hospital bed Mattress, Over-bed table, Beside cabinet</p>
	two functions hospital patient bed	AC-MB017	<p>Specification: L2150*W950*H500mm</p> <ol style="list-style-type: none"> 1. Detachable ABS head/foot board 3. Integrated stamping frame and mattress base are made of steel epoxy 4. I.V. Rod & Drainage hook 5. 5"Castor with two brakes 6. Optional accessory: Mattress & over-bed table & bedside cabinet <p>Back-rest lifting 0-75°(±2%) Knee-rest lifting 0-45°(±2%) Optional Accessory: second hand hospital bed second hand hospital bed second hand hospital bed Mattress, Over-bed table, Beside cabinet</p>
	Single function hospital bed cost	AC-MB018	<p>Specification: L2150*W950*H500mm</p> <p>Standard Accessory Detachable ABS head/foot board Collapsible aluminum alloy guardrails</p> <p>Integrated stamping frame and mattress base are made of steel epoxy I V Pole& Drainage hook with PU cover Heavy duty caster with brakes 5 inch Manual crank system 1set Optional accessory: Mattress over-bed table bedside cabinet</p>
	Single function hospital bed cost	AC-MB019	<p>Related parameters: Overall dimension: L2130*W980*H490mm</p> <p>Standard configuration: 1 set manual crank systems: 1 pair plastic head/foot-end boards; 1 pair Aluminium siderails 4pcs urine hooks; 1pc S.S I.V pole;</p>

	Flat bed	AC-MB020	<p>Related parameters: Overall dimension: L2130*W980*H490mm Standard configuration: 1 pair plastic head/foot-end boards; 4pcs urine hooks; 1pc S.S I.V pole;</p>
	Single function hospital bed cost	AC-MB021	<p>Related parameters: Overall dimension: L2130*W980*H490mm Standard configuration: 1 set manual crank systems; 1 pair plastic head/foot-end boards; 4pcs urine hooks; 1pc S.S I.V pole;</p>

	Orthopedics Traction Bed	AC-MB022	<p>Specification: L2150*W950*H500mm Standard Accessory: Detachable ABS head/foot board Collapsible aluminum alloy guardrails Integrated stamping frame and mattress base are made of steel epoxy Traction frame by stainless steel I V Pole Drainage hook with PU cover Heavy duty caster with brakes 5" Manual crank system 2set Optional accessory: Mattress over-bed table bedside cabinet</p>
---	-----------------------------	----------	---

	Orthopedics Traction Bed	AC-MB023	<p>Specification L2150*950*H500mm Standard Accessory Detachable ABS head/foot board Collapsible aluminum alloy guardrails Integrated stamping frame and mattress base are made of steel epoxy OT raction frame by stainless steel I V Pole Drainage hook with PU cover Heavy duty caster with brakes 5 Manual crank system 2set Optional accessory Mattress over-bed table bedside cabinet</p>
	Four cranks manual Orthopedics Traction Bed	AC-MB024	<p>Related parameters: Overall dimension: L2200*W980*H490-720mm Back rest tilting: 80°±5°, knee board tilting: 40°±5° Hi-lo adjustment between:490-720mm Standard configuration: 4 sets manual crank systems; 4pcs luxurious castors; 1 pair plastic head/foot-end boards; 1 pair Aluminium siderails 4pcs urine hooks; 1pc S.S I.V pole; 1 set orthopedics traction frame; Options: Medical mattress; dining board; bedside cabinet; patient file panel; sundries basket</p>
	Four cranks manual Orthopedics Traction Bed	AC-MB025	<p>Related parameters: Overall dimension: L2200*W980*H490-720mm Back rest tilting: 80°±5°, knee board tilting: 40°±5° Hi-lo adjustment between:490-720mm Standard configuration: 4 sets manual crank systems; central-controlled brake system,4pcs luxurious cnetral locking castors; 1 pair plastic head/foot-end boards; 1 pair Aluminium siderails 4pcs urine hooks; 1pc S.S I.V pole; 1 set orthopedics traction frame; Options: Medical mattress; dining board; bedside cabinet; patient file panel; sundries basket</p>
	Three cranks manual Orthopedics Traction Bed	AC-MB026	<p>Related parameters: Overall dimension: L2200*W980*H490mm Back rest tilting: 80°±5°, knee board tilting: 40°±5° Standard configuration: 3 sets manual crank systems; 4pcs luxurious castors; 1 pair plastic head/foot-end boards; 1 pair Aluminium siderails 4pcs urine hooks; 1pc S.S I.V pole; 1 set orthopedics traction frame; Options: Medical mattress; dining board; bedside cabinet; patient file panel; sundries basket</p>