














ZHANGJIAGANG ANNECY MACHINERY CO., LTD.

MEDICAL WORKSTATION TROLLEY CATALOGUE


Images	Description	Model No	Details	SIZE
	Workstation trolley	AC-WNT001	Tabletop is made of ABS material by injection molding techniques; Height of tabletop adjustable; With space for mouse,keyboard; Two drawers; Fan-shaped chassis,move more steady	480*520*830-1030MM ,Packing:650*690*1082mm
	Workstation trolley	AC-WNT002	Tabletop is made of ABS material by injection molding techniques; Height of tabletop adjustable; With space for mouse,keyboard; Two drawers; Fan-shaped chassis,move more steady	480*520*830-1030MM ,Packing:650*690*1082mm
	Workstation trolley	AC-WNT003	Tabletop is made of ABS material by injection molding techniques; Height of tabletop adjustable; With space for mouse,keyboard; Two drawers; Fan-shaped chassis,move more steady	480*520*830-1030MM ,Packing:650*690*1082mm
	Workstation trolley	AC-WNT004	Tabletop is made of ABS material by injection molding techniques; Height of tabletop adjustable; With space for mouse,keyboard; Fan-shaped chassis,move more steady	480*520*830-1030MM ,Packing:650*690*1082mm

	Workstation trolley	AC-WNT005	Tabletop is made of ABS material by injection molding techniques; Height of tabletop adjustable; With space for mouse,keyboard; Shelf to keep printer; Fan-shaped chassis,move more steady	480*520*830-1030MM ,Packing:650*690*1082mm
	Workstation trolley	AC-WNT006	Tabletop is made of ABS material by injection molding techniques; Height of tabletop adjustable; With space for mouse,keyboard; Fan-shaped chassis,move more steady	480*520*830-1030MM ,Packing:650*690*1082mm
	Workstation trolley	AC-WNT007	Tabletop is made of ABS material by injection molding techniques; Height of tabletop adjustable; One bucket; Fan-shaped chassis,move more steady	480*520*830-1080MM
	Workstation trolley	AC-WNT008	Tabletop is made of ABS material by injection molding techniques; Fan-shaped chassis,move more steady	480*520*920mm
	Workstation trolley	AC-WNT009	Tabletop is made of ABS material by injection molding techniques; Two buckets; Fan-shaped chassis,move more steady	480*520*920mm


	Workstation trolley	AC-WNT010	Tabletop is made of ABS material by injection molding techniques; One shelf ; Fan-shaped chassis,move more steady	480*520*920mm
	Workstation trolley	AC-WNT011	Tabletop is made of ABS material by injection molding techniques; One shelf ; One drawer; Fan-shaped chassis,move more steady	480*520*920mm
	Workstation trolley	AC-WNT012	Tabletop is made of ABS material by injection molding techniques; One shelf ; Two drawers; Fan-shaped chassis,move more steady	480*520*920mm
	Workstation trolley	AC-WNT013	Tabletop is made of ABS material by injection molding techniques; One shelf ; Fan-shaped chassis,move more steady	480*520*920mm
	Workstation trolley	AC-WNT014	Tabletop is made of ABS material by injection molding techniques; Two shelves ; One drawer; Fan-shaped chassis,move more steady	480*520*920mm
	Workstation trolley	AC-WNT015	Tabletop is made of ABS material by injection molding techniques; Three shelves; Fan-shaped chassis,move more steady	480*520*920mm


	Workstation trolley	AC-WNT016	<p>Tabletop is made of ABS material by injection molding techniques; Built-in host computer,low power loss,stabilized D.C source; Lithium battery;12-18.5V,long battery life and safe; With space for mouse,keyboard; Three drawers; It is safe to use and to manage to fix the battery upside; Display holder:applied for ≤22" display,proposal load ≤10kg</p>	520*510*1020mm
	Workstation trolley	AC-WNT017	<p>Tabletop is made of ABS material by injection molding techniques; Five drawers; Display holder:applied for ≤22" display,proposal load ≤10kg</p>	520*510*1020mm
	Workstation trolley	AC-WNT018	<p>Tabletop is made of ABS material by injection molding techniques; Built-in host computer,low power loss,stabilized D.C source; Lithium battery;12-18.5V,long battery life and safe; With space for mouse,keyboard; Five drawers; It is safe to use and to manage to fix the battery upside; Display holder:applied for ≤22" display,proposal load ≤10kg</p>	520*510*1020mm
	Workstation trolley	AC-WNT019	<p>Tabletop is made of ABS material by injection molding techniques; With space for mouse,keyboard; Five drawers; It is safe to use and to manage to fix the battery upside; Display holder:applied for ≤22" display,proposal load ≤10kg</p>	520*510*1020mm
	Workstation trolley	AC-WNT020	<p>Tabletop is made of ABS material by injection molding techniques; Three drawers; One shelf</p>	520*510*1020mm
	Workstation trolley	AC-WNT021	<p>Tabletop is made of ABS material by injection molding techniques; Built-in host computer,low power loss,stabilized D.C source; Lithium battery;12-18.5V,long battery life and safe; With space for mouse,keyboard; Two drawers; It is safe to use and to manage to fix the battery upside; Display holder:applied for ≤22" display,proposal load ≤10kg</p>	520*510*1020mm


	Workstation trolley	AC-WNT022	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>With space for mouse,keyboard;</p> <p>Three drawers;</p> <p>It is safe to use and to manage to fix the battery upside;</p> <p>Display holder:applied for ≤22" display,proposal load ≤10kg</p>	520*510*1020mm
	Workstation trolley	AC-WNT023	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Four drawers;</p> <p>One shelf</p> <p>Display holder:applied for ≤22" display,proposal load ≤10kg</p> 	780*475*1000mm
	Workstation trolley	AC-WNT024	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Six drawers;</p> <p>One shelf</p> <p>Display holder:applied for ≤22" display,proposal load ≤10kg</p>	780*475*1000mm
	Workstation trolley	AC-WNT025	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Four drawers;</p> <p>Display holder:applied for ≤22" display,proposal load ≤10kg</p>	628*475*930mm
	Workstation trolley	AC-WNT026	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Four drawers;</p> <p>Display holder:applied for ≤22" display,proposal load ≤10kg</p> 	628*475*930mm
	Workstation trolley	AC-WNT027	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Five drawers;</p> <p>Display holder:applied for ≤22" display,proposal load ≤10kg</p>  	626*475*930mm


	Workstation trolley	AC-WNT028	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Five drawers;</p> <p>Display holder:applied for $\leq 22^{\circ}$ display,proposal load $\leq 10\text{kg}$</p>	760*475*930mm
---	---------------------	-----------	--	---------------





	Workstation trolley	AC-WNT029	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Five drawers;</p> <p>Display holder:applied for $\leq 22^{\circ}$ display,proposal load $\leq 10\text{kg}$</p>	625*475*930mm
--	---------------------	-----------	--	---------------

	Workstation trolley	AC-WNT030	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Two drawers; Two doors;</p> <p>Display holder:applied for $\leq 22^{\circ}$ display,proposal load $\leq 10\text{kg}$</p>	625*475*930mm
---	---------------------	-----------	--	---------------

	Workstation trolley	AC-WNT031	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Four drawers; One shelf</p> <p>Display holder:applied for $\leq 22^{\circ}$ display,proposal load $\leq 10\text{kg}$</p>	750*475*930mm
---	---------------------	-----------	--	---------------

	Workstation trolley	AC-WNT032	<p>Tabletop is made of ABS material by injection molding techniques;</p> <p>Five drawers;</p> <p>Display holder:applied for $\leq 22^{\circ}$ display,proposal load $\leq 10\text{kg}$</p>	625*475*930mm
---	---------------------	-----------	--	---------------

	Workstation trolley	AC-WNT033	Tabletop is made of ABS material by injection molding techniques; Five drawers; Display holder:applied for ≤22" display,proposal load ≤10kg	625*475*930mm
	Workstation trolley	AC-WNT034	Tabletop is made of ABS material by injection molding techniques; Five drawers; Display holder:applied for ≤22" display,proposal load ≤10kg	750*475*930mm