






**HANGJIAGANG ANNECY MACHINERY CO., LTD.**

**Medicine trolley Catalogue & Quotation**


Images	Model No	Details	Size
	AC-MT001	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers	850*520*960mm
	AC-MT002	2.3*3 segregation board in every drawer 3.double trash can 4.Side tumbling type material box	750*475*930mm
	AC-MT003	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers	750*475*930mm
	AC-MT004	2.3*3 segregation board in every drawer 3.double trash can 4.Side tumbling type material box	850*520*960mm
	AC-MT005	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.50 person double face medicine supply 2.3*3 segregation board in every drawer 3.Double trash can 4.Side tumbling type material box	850*520*1010mm
	AC-MT006	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.50 person double face medicine supply 2.3*3 segregation board in every drawer 3.Double trash can 4.Side tumbling type material box	750*475*930mm
         			
	AC-MT007	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Six drawers	850*520*960mm
	AC-MT008	2.3*3 segregation board in every drawer 3.double trash can 4.Side tumbling type material box	750*475*930mm
	AC-MT009	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1 Four drawers	673*480*960mm

	AC-MT010	1.1. Storage basket 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	625*475*930mm
	AC-MT011	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Six drawers 3*3 segregation board in every drawer 2.Storage box 3.One big trash can one small dust can	630*470*930mm
	AC-MT012	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	673*480*960mm



	AC-MT013	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	688*480*960mm
	AC-MT014		630*470*930mm
	AC-MT015	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	850*520*960mm
	AC-MT016		720*470*930mm



	AC-MT017	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration:	850*520*960mm
---	----------	--	---------------

	AC-MT018	Configuration: 1.Six drawers 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	720*470*930mm
			
	AC-MT019	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage box 3.3*3 segregation board in every drawer 4.Double trash can	850*520*960mm
	AC-MT020		720*470*930mm
			
	AC-MT021	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage box 3.3*3 segregation board in every drawer 4.Double trash can	850*520*960mm
	AC-MT022		720*470*930mm
			
	AC-MT023	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.One small drawer and one middle drawer 2.3*3 segregation board in every drawer 3.double trash can	750*475*930mm
	AC-MT024	4.Side tumbling type material box 5.Cabinet with two doors	850*520*960mm

	AC-MT025	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	750*475*930mm
			
	AC-MT026	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	750*475*930mm
	AC-MT027		850*520*930mm
			
 	AC-MT028	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Six drawers 2.3*3 segregation board in every drawer 3.double trash can 4.Side tumbling type material box	750*475*930mm
	AC-MT029		850*520*930mm
	AC-MT030	Mainly made of plastic, aluminum and steel ABS top board ,easy to clean Configuration: 1.Four drawers 2.Storage basket 3.3*3 segregation board in every drawer 4.Double trash can	690*482*930mm
	AC-MT031		623*475*930mm

