






























# ZHANGJIAGANG ANNECY MACHINERY CO., LTD.


## CATALOGUE & QUOTATION










Images	Description	Model No	Details	Packing
	Shadowless Operating lamp	AC-OL001	Illuminance: $\geq 140,000\text{Lux}$ $\geq 140,000\text{Lux}$ Colour temperature: $4300\pm 200\text{K}$ Colour reduction index(Ra):93 Illuminance depth: $\geq 1300\text{mm}$ Total irradiance: $426\text{W/m}^2$ The ratio between irradiance and illuminance : $3.3\text{mW/m}^2\text{lx}$ Size of light field: $160 - 280\text{mm}$ LED bulb: $37\text{V}/50\text{W}$ Service life of illuminant: $30,000\text{h}$ Power supply voltage: $\text{AC}110-240\text{v}, 50/60\text{Hz}$ Brightness adjustment: Automatic 8-stage continuous light adjustment Lowest height of installation: $2900\text{mm}$	2 piece in a carton: $91^*92^*37$ (cm) GW: $19\text{kg}$ NW: $12\text{kg}$ 1 piece in a wooden case: $122^*43^*72$ (cm) GW: $65\text{kg}$ NW: $50\text{kg}$
	Shadowless Operating lamp	AC-OL002	Illuminance $\geq 140,000\text{Lux}$ Colour reduction Colour temperature: $4300\pm 200\text{K}$ index(Ra):93 Illuminance depth: $\geq 1300\text{mm}$ Total irradiance: $426\text{W/m}^2$ The ratio between irradiance and illuminance : $3.3\text{mW/m}^2\text{lx}$ Size of light field: $160 - 280\text{mm}$ LED bulb: $37\text{V}/50\text{W}$ Service life of illuminant: $30,000\text{h}$ Power supply voltage: $\text{AC}110-240\text{v}, 50/60\text{Hz}$ Brightness adjustment: Automatic 8-stage continuous light adjustment Lowest height of installation: $2900\text{mm}$	1 piece in a carton: $91^*92^*37$ (cm) GW: $19\text{kg}$ NW: $12\text{kg}$ 1 piece in a wooden case: $122^*43^*72$ (cm) GW: $65\text{kg}$ NW: $50\text{kg}$
	LED Shadowless Operating lamp	AC-OL003	Illuminance: $\geq 160,000\text{Lux}$ $\geq 160,000\text{Lux}$ Colour temperature: $4500\pm 200\text{K}$ Colour reduction index(Ra):96 Illuminance depth: $1150\text{mm}/1150\text{mm}$ Total irradiance: $510\text{W/m}^2$ Size of light field: $180-260\text{mm}$ lamp head diameter: $760\text{mm}/760\text{mm}$ Rated power of bulb: $3.2\text{V}/1\text{W}$ Service life of illuminant: $30,000\text{h}$ Power supply voltage: $\text{AC}110-240\text{v}, 50/60\text{Hz}$ Lowest height of installation: $2900\text{mm}$	2 piece in a carton: $113^*109^*45$ (cm) GW: $36\text{kg}$ NW: $27\text{kg}$ 1 piece in a wooden case: $122^*42^*85$ (cm) GW: $81\text{kg}$ NW: $70\text{kg}$
	LED Shadowless Operating lamp	AC-OL004	Illuminance: $\geq 160,000\text{Lux}$ Colour reduction Colour temperature: $4500\pm 200\text{K}$ index(Ra):96 Illuminance depth: $1150\text{mm}/1150\text{mm}$ Total irradiance: $510\text{W/m}^2$ Size of light field: $180-260\text{mm}$ lamp head diameter: $760\text{mm}/760\text{mm}$ Rated power of bulb: $3.2\text{V}/1\text{W}$ Service life of illuminant: $30,000\text{h}$ Power supply voltage: $\text{AC}110-240\text{v}, 50/60\text{Hz}$ Lowest height of installation: $2900\text{mm}$	1 piece in a carton: $113^*109^*45$ (cm) GW: $36\text{kg}$ NW: $27\text{kg}$ 1 piece in a wooden case: $122^*42^*85$ (cm) GW: $81\text{kg}$ NW: $70\text{kg}$
	EM-SL01-LED5+5 - TV (LED) Shadowless Operating lamp with camera	AC-OL005	Illuminance: $\geq 160,000\text{Lux}$ $\geq 160,000\text{Lux}$ Colour temperature: $4500\pm 200\text{K}$ Colour reduction index(Ra):96 Illuminance depth: $1150\text{mm}/1150\text{mm}$ Total irradiance: $510\text{W/m}^2$ Size of light field: $200 - 350\text{mm}$ lamp head diameter: $860\text{mm}/860\text{mm}$ Rated power of bulb: $3.2\text{V}/1\text{W}$ Service life of illuminant: $30,000\text{h}$ Power supply voltage: $\text{AC}110-240\text{v}, 50/60\text{Hz}$ Lowest height of installation: $2900\text{mm}$	2 piece in a carton: $113^*109^*45$ (cm) GW: $36\text{kg}$ NW: $27\text{kg}$ 1 piece in a wooden case: $122^*42^*85$ (cm) GW: $81\text{kg}$ NW: $70\text{kg}$
	EM-SL01-LED5+5 (LED) Shadowless Operating lamp	AC-OL006	Illuminance: $\geq 160,000\text{Lux}$ $\geq 160,000\text{Lux}$ Colour temperature: $4500\pm 200\text{K}$ Colour reduction index(Ra):96 Illuminance depth: $1150\text{mm}/1150\text{mm}$ Total irradiance: $510\text{W/m}^2$ Size of light field: $200 - 350\text{mm}$ lamp head diameter: $860\text{mm}/860\text{mm}$ Rated power of bulb: $3.2\text{V}/1\text{W}$ Service life of illuminant: $30,000\text{h}$ Power supply voltage: $\text{AC}110-240\text{v}, 50/60\text{Hz}$ Lowest height of installation: $3100\text{mm}$ Camera data:	2 piece in a carton: $113^*109^*45$ (cm) GW: $36\text{kg}$ NW: $27\text{kg}$ 1 piece in a wooden case: $122^*42^*85$ (cm) GW: $81\text{kg}$ NW: $70\text{kg}$ 1 piece in a wooden
	EM-SL02-LED3+5 (LED) Shadowless Operating lamp	AC-OL007	Illuminance: $\text{LED}3 \geq 100,000\text{Lux}$ LED5: $\geq 160,000\text{Lux}$ Colour temperature: $4500\pm 200\text{K}$ Colour reduction index(Ra):96 Illuminance depth: $1020\text{mm}/1150\text{mm}$ Rated power of bulb: $3.2\text{V}/1\text{W}$ Total irradiance: $\text{LED}3: 380\text{W/m}^2$ Size of light field: $200 - 350\text{mm}$ LED5: $510\text{W/m}^2$ Service life of illuminant: $30,000\text{h}$ Power supply voltage: $\text{AC}110-240\text{v}, 50/60\text{Hz}$ Lowest height of installation: $2900\text{mm}$	1 piece in a carton: $113^*109^*45$ (cm) GW: $28\text{kg}$ NW: $20\text{kg}$ 1 piece in a wooden case: $104^*97^*44$ (cm) GW: $36\text{kg}$ NW: $27\text{kg}$ 1 piece in a
	EM-SL03-LED3+3 (LED) Shadowless Operating lamp	AC-OL008	Illuminance: $\geq 100,000\text{Lux}$ $\geq 100,000\text{Lux}$ Colour temperature: $4500\pm 200\text{K}$ Colour reduction index:96 Illuminance depth: $1020\text{mm}/1020\text{mm}$ Rated power of bulb: $3.2\text{V}/1\text{W}$ Total irradiance: $380\text{W/m}^2$ Size of light field: $200 - 350\text{mm}$ Service life of illuminant: $30,000\text{h}$ Lowest height of installation: $2900\text{mm}$	2 piece in a carton: $104^*97^*44$ (cm) GW: $36\text{kg}$ NW: $27\text{kg}$ 1 piece in a wooden case: $122^*42^*85$ (cm) GW: $81\text{kg}$ NW: $70\text{kg}$










	EM-SL04-LED5 (LED) Shadowless Operating lamp	AC-OL009	<p>           Illuminance: <math>\geq 160,000</math> Lux            Colour temperature: <math>4500 \pm 200</math> K            Colour reduction index (Ra): 96            power of bulb: 3.2V/1W            Total irradiance: <math>510</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            Lowest height of installation: 2900 mm         </p> <p>           Illuminance depth: 1150 mm            Rated         </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz         </p>	<p>           1 piece in a carton:            113*109*45 (cm)            GW: 36kg            NW: 27kg            1 piece in a wooden case:            110*42*72 (cm)            GW: 49kg            NW: 38kg         </p>
	EM-SL05-LED3W (LED) Shadowless Operating lamp	AC-OL010	<p>           Illuminance: <math>\geq 100,000</math> Lux            Colour temperature: <math>4500 \pm 200</math> K            Colour reduction index (Ra): 96            power of bulb: 3.2V/1W            Total irradiance: <math>380</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            Lowest height of installation: 2900 mm         </p> <p>           Illuminance depth: 1020 mm            Rated         </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz         </p>	<p>           1 piece in a carton: 104*97*44 (cm)            GW: 36kg            NW: 27kg            1 piece in a wooden case:            110*42*72 (cm)            GW: 49kg            NW: 38kg         </p>
	EM-SL05-LED3 (LED) Shadowless Operating lamp	AC-OL011	<p>           Illuminance: <math>\geq 100,000</math> Lux            Colour temperature: <math>4500 \pm 200</math> K            Colour reduction index (Ra): 96            power of bulb: 3.2V/1W            Total irradiance: <math>380</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            Lowest height of installation: 2900 mm         </p> <p>           Illuminance depth: 1020 mm            Rated         </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz         </p>	<p>           1 piece in a carton: 104*97*44 (cm)            GW: 36kg            NW: 27kg            1 piece in a wooden case:            110*42*72 (cm)            GW: 49kg            NW: 38kg         </p>
	EM-SL06-LED3S(on stand) (LED) Shadowless Operating lamp	AC-OL012	<p>           Illuminance: <math>\geq 100,000</math> Lux            Colour temperature: <math>4500 \pm 200</math> K            Colour reduction index (Ra): 96            Power supply voltage: AC110-240v, 50/60 Hz            Total irradiance: <math>380</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h         </p> <p>           Illuminance depth: 1020 mm            Rated power of bulb: 3.2V/1W         </p>	<p>           1 piece in a wooden case:            106*98*51 (cm)            GW: 55kg            NW: 42kg            1 piece in a wooden case:            86*51*24 (cm)            GW: 36kg            NW: 30kg         </p>
	EM-SL06-LED3E(on stand) (Emergency) (LED) Shadowless Operating lamp	AC-OL013	<p>           Illuminance: <math>\geq 100,000</math> Lux            Colour temperature: <math>4500 \pm 200</math> K            Colour reduction index (Ra): 96            Rated power of bulb: 3.2V/1W            Total irradiance: <math>380</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            240v, 50/60 Hz         </p> <p>           Illuminance depth: 1020 mm         </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz            Emergency time 2h         </p>	<p>           1 piece in a wooden case:            106*98*51 (cm)            GW: 55kg            NW: 42kg            1 piece in a wooden case:            76*53*38 (cm)            GW: 44kg            NW: 37kg         </p>
	EM-SL07-LED5+5 (Adjust color temperature) (LED) Shadowless Operating lamp	AC-OL014	<p>           Illuminance: <math>\geq 140,000</math> Lux            Colour temperature: <math>3500 \pm 200</math> K/<math>4000 \pm 200</math> K            Colour reduction index (Ra): 96            depth: 1150 mm/1150 mm            Total irradiance: <math>510</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            Lowest height of installation: 2900 mm         </p> <p>           Illuminance            Rated power of bulb: 3.2V/1W         </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz         </p>	<p>           2 piece in a carton:            113*109*45 (cm)            GW: 36kg            NW: 27kg            1 piece in a wooden case:            122*42*85 (cm)            GW: 81kg            NW: 70kg         </p>
	EM-SL08-LED3+5 (Adjust color temperature) (LED) Shadowless Operating lamp	AC-OL015	<p>           Illuminance: LED3: <math>\geq 100,000</math> Lux            Colour temperature: <math>3500 \pm 200</math> K/<math>4000 \pm 200</math> K            Colour reduction index (Ra): 96            Total irradiance: LED3: <math>380</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            Lowest height of installation: 2900 mm         </p> <p>           LED5: <math>\geq 140,000</math> Lux            4500 <math>\pm</math> </p> <p>           LED5: <math>510</math> W/m<sup>2</sup> </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz         </p>	<p>           1 piece in a carton:            113*109*45 (cm)            GW: 28kg            NW: 20kg            1 piece in a carton:            104*97*44 (cm)            GW: 36kg            NW: 27kg            1 piece in a         </p>
	EM-SL09-LED3+3 (Adjust color temperature) (LED) Shadowless Operating lamp	AC-OL016	<p>           Illuminance: <math>\geq 100,000</math> Lux            Colour temperature: <math>3500 \pm 200</math> K/<math>4000 \pm 200</math> K            Colour reduction index (Ra): 96            Total irradiance: <math>380</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            Lowest height of installation: 2900 mm         </p> <p>           Illuminance depth: 1020 mm/1020 mm         </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz         </p>	<p>           2 piece in a carton: 104*97*44 (cm)            GW: 36kg            NW: 27kg            1 piece in a wooden case:            122*42*85 (cm)            GW: 81kg            NW: 70kg         </p>
	EM-SL10-LED5 (Adjust color temperature) (LED) Shadowless Operating lamp	AC-OL017	<p>           Illuminance: <math>\geq 140,000</math> Lux            Colour temperature: <math>3500 \pm 200</math> K/<math>4000 \pm 200</math> K/<math>4500 \pm 200</math> K/<math>5000 \pm 200</math> K            Colour reduction index (Ra): 96            Total irradiance: <math>510</math> W/m<sup>2</sup>            Size of light field: 200 – 350 mm            Service life of illuminant: 30,000 h            Lowest height of installation: 2900 mm         </p> <p>           Illuminance depth: 1150 mm         </p> <p>           Power supply voltage: AC110-240v, 50/60 Hz         </p>	<p>           1 piece in a carton:            113*109*45 (cm)            GW: 36kg            NW: 27kg            1 piece in a wooden case:            110*42*72 (cm)            GW: 49kg            NW: 38kg         </p>










	EM-SL11-LED3 (Adjust color temperature) (LED)Shadowless Operating lamp	AC-OL018	<p>Illuminance:≥100,000Lux Colour temperature:3500±200K/4000±200K 4500±200K/5000±200K Colour reduction index(Ra):96 Total irradiance:380W/m<sup>2</sup> Size of light field:200~350mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>Illuminance depth:1020mm</p> <p>Power supply voltage:AC110-240v,50/60Hz</p>	<p>1 piece in a carton: 104*97*44 (cm) GW:36kg NW:27kg 1 piece in a wooden case:110*42*72 (cm) GW:49kg NW:38kg</p>
	EM-SL12-LED3S(on stand) (Adjust color temperature) (LED)Shadowless Operating lamp	AC-OL019	<p>Illuminance:≥100,000Lux Colour temperature:3500±200K/4000±200K 200K 4500±200K/5000±200K Colour reduction index(Ra):96 Total irradiance:380W/m<sup>2</sup> Size of light field:200~350mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>Illuminance depth:1020mm</p> <p>Power supply voltage:AC110-240v,50/60Hz</p>	<p>1 piece in a wooden case: 106*98*51 (cm) GW:55kg NW:42kg 1 piece in a wooden case: 86*51*24 (cm) GW:36kg NW:30kg</p>
	EM-SL12-LED3E(on stand)(Emergency) (Adjust color temperature) (LED)Shadowless Operating lamp	AC-OL020	<p>Illuminance:≥100,000Lux Colour temperature:3500±200K/4000±200K 200K 4500±200K/5000±200K Colour reduction index(Ra):96 Total irradiance:380W/m<sup>2</sup> Size of light field:200~350mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>Illuminance depth:1020mm</p> <p>Power supply voltage:AC110-240v,50/60Hz Emergency time 2h</p>	<p>1 piece in a wooden case: 106*98*51 (cm) GW:55kg NW:42kg 1 piece in a wooden case: 76*53*38 (cm) GW:44kg NW:37kg</p>
	LED Shadowless Operating lamp	AC-OL021	<p>Illuminance:≥130,000Lux Colour temperature:4800±200K Colour reduction index(Ra):96 Total irradiance:480W/m<sup>2</sup> Size of light field:120~240mm diameter:550mm/550mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>≥130,000Lux Illuminance depth:1050mm</p> <p>lamp head Rated power of bulb:3.2V/1W Power supply voltage:AC110-240v,50/60Hz</p>	<p>2 piece in a carton: 91*92*37 (cm) GW:21kg NW:14kg 1 piece in a wooden case: 122*43*72 (cm) GW:65kg NW:50kg</p>
	LED Shadowless Operating lamp	AC-OL022	<p>Illuminance:≥120,000Lux Colour temperature:4800±200K Colour reduction index(Ra):96 Total irradiance:420W/m<sup>2</sup> Size of light field:120~240mm diameter:550mm/550mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>≥120,000Lux Illuminance depth:1050mm</p> <p>lamp head Rated power of bulb:3.2V/1W Power supply voltage:AC110-240v,50/60Hz</p>	<p>2 piece in a carton: 80*81*36 (cm) GW:19kg NW:13kg 1 piece in a wooden case: 122*43*72 (cm) GW:65kg NW:50kg</p>
	LED Shadowless Operating lamp	AC-OL023	<p>Illuminance:≥100,000Lux Colour temperature:4800±200K Colour reduction index(Ra):96 Total irradiance:310W/m<sup>2</sup> Size of light field:120~240mm diameter:500mm/500mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>≥100,000Lux Illuminance depth:1050mm</p> <p>lamp head Rated power of bulb:3.2V/1W Power supply voltage:AC110-240v,50/60Hz</p>	<p>2 piece in a carton: 80*81*36 (cm) GW:16kg NW:10kg 1 piece in a wooden case: 122*43*72 (cm) GW:65kg NW:50kg</p>
	LED Shadowless Operating lamp	AC-OL024	<p>Illuminance:LED3:≥100,000Lux Colour temperature:4800±200K Colour reduction index(Ra):96 Total irradiance:LED3:310W/m<sup>2</sup> Size of light field:120~240mm diameter:500mm/550mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>LED4:≥120,000Lux Illuminance depth:1050mm LED4:420W/m<sup>2</sup> lamp head Rated power of bulb:3.2V/1W Power supply voltage:AC110-240v,50/60Hz</p>	<p>1 piece in a carton: 80*81*36 (cm) GW:16kg NW:10kg 1 piece in a carton: 80*81*36 (cm) GW:19kg NW:13kg 1 piece in a wooden case:</p>
	LED Shadowless Operating lamp	AC-OL025	<p>Illuminance:LED3:≥100,000Lux Colour temperature:4800±200K Colour reduction index(Ra):96 Total irradiance:LED3:310W/m<sup>2</sup> Size of light field:120~240mm diameter:500mm/550mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>LED5:≥130,000Lux Illuminance depth:1050mm LED5:480W/m<sup>2</sup> lamp head Rated power of bulb:3.2V/1W Power supply voltage:AC110-240v,50/60Hz</p>	<p>1 piece in a carton: 80*81*36 (cm) GW:16kg NW:10kg 1 piece in a carton: 91*92*37 (cm) GW:21kg NW:14kg 1 piece in a wooden case:</p>
	LED Shadowless Operating lamp	AC-OL026	<p>Illuminance:LED4:≥120,000Lux Colour temperature:4800±200K Colour reduction index(Ra):96 Total irradiance:LED4:420W/m<sup>2</sup> Size of light field:120~240mm diameter:550mm/550mm Service life of illuminant:30,000h Lowest height of installation:2900mm</p> <p>LED5:≥130,000Lux Illuminance depth:1050mm LED5:480W/m<sup>2</sup> lamp head Rated power of bulb:3.2V/1W Power supply voltage:AC110-240v,50/60Hz</p>	<p>1 piece in a carton: 80*81*36 (cm) GW:19kg NW:13kg 1 piece in a carton: 91*92*37 (cm) GW:21kg NW:14kg 1 piece in a wooden case:</p>

	LED Shadowless Operating lamp	AC-OL027	<p>           Illuminance:≥130,000Lux            Colour temperature:4800±200K            Colour reduction index(Ra):96            Total irradiance:480W/m<sup>2</sup>            Size of light field:120~240mm            Rated power of bulb:3.2V/1W            Service life of illuminant:30,000h            240v,50/60Hz            Lowest height of installation:2900mm         </p>	<p>           Illuminance depth:1050mm            lamp head diameter:550mm            Power supply voltage:AC110-240v,50/60Hz         </p>	<p>           1 piece in a carton: 91*92*37 (cm)            GW:21kg            NW:14kg            1 piece in a wooden case: 108*43*60 (cm)            GW:65kg            NW:50kg         </p>
	LED Shadowless Operating lamp	AC-OL028	<p>           Illuminance:≥100,000Lux            Colour temperature:4800±200K            Colour reduction index(Ra):96            Total irradiance:420W/m<sup>2</sup>            Size of light field:120~240mm            Rated power of bulb:3.2V/1W            Service life of illuminant:30,000h            240v,50/60Hz            Lowest height of installation:2900mm         </p>	<p>           Illuminance depth:1050mm            lamp head diameter:550mm            Power supply voltage:AC110-240v,50/60Hz         </p>	<p>           1 piece in a carton: 80*81*36 (cm)            GW:19kg            NW:13kg            1 piece in a wooden case: 108*43*60 (cm)            GW:65kg            NW:50kg         </p>
	LED Shadowless Operating lamp	AC-OL029	<p>           Illuminance:≥80,000Lux            Colour temperature:4800±200K            Colour reduction index(Ra):96            Total irradiance:310W/m<sup>2</sup>            Size of light field:120~240mm            Rated power of bulb:3.2V/1W            Service life of illuminant:30,000h            240v,50/60Hz            Lowest height of installation:2900mm         </p>	<p>           Illuminance depth:1050mm            lamp head diameter:500mm            Power supply voltage:AC110-240v,50/60Hz         </p>	<p>           1 piece in a carton: 80*81*36 (cm)            GW:16kg            NW:10kg            1 piece in a wooden case: 108*43*60 (cm)            GW:65kg            NW:50kg         </p>
	LED Shadowless Operating lamp	AC-OL030	<p>           Illuminance:≥130,000Lux            Colour temperature:4800±200K            Colour reduction index(Ra):96            Total irradiance:480W/m<sup>2</sup>            Size of light field:120~240mm            diameter:550mm            Rated power of bulb:3.2V/1W            Service life of illuminant:30,000h            240v,50/60Hz         </p>	<p>           Illuminance depth:1050mm            lamp p head            Power supply voltage:AC110-240v,50/60Hz         </p>	<p>           1 piece in a wooden case: 95*92*44 (cm)            GW:41kg            NW:27kg            1 piece in a wooden case: 86*51*24 (cm)            GW:36kg            NW:30kg         </p>
	LED Shadowless Operating lamp	AC-OL031	<p>           Illuminance:≥120,000Lux            Colour temperature:4800±200K            Colour reduction index(Ra):96            Total irradiance:420W/m<sup>2</sup>            Size of light field:120~240mm            diameter:550mm            Rated power of bulb:3.2V/1W            Service life of illuminant:30,000h            240v,50/60Hz         </p>	<p>           Illuminance depth:1050mm            lamp p head            Power supply voltage:AC110-240v,50/60Hz         </p>	<p>           1 piece in a wooden case: 95*92*44 (cm)            GW:40kg            NW:26kg            1 piece in a wooden case: 86*51*24 (cm)            GW:36kg            NW:30kg         </p>
	LED Shadowless Operating lamp	AC-OL032	<p>           Illuminance:≥100,000Lux            Colour temperature:4800±200K            Colour reduction index(Ra):96            Total irradiance:310W/m<sup>2</sup>            Size of light field:120~240mm            Rated power of bulb:3.2V/1W            Service life of illuminant:30,000h            240v,50/60Hz         </p>	<p>           Illuminance depth:1050mm            lamp head diameter:500mm            Power supply voltage:AC110-240v,50/60Hz            Emergency time: 2.5 hours         </p>	<p>           1 piece in a wooden case: 95*92*44 (cm)            GW:37kg            NW:23kg            1 piece in a wooden case: 86*51*24 (cm)            GW:36kg            NW:30kg         </p>
	on stand LED Shadowless Operating lamp	AC-OL033	<p>           Illuminance:≥100,000Lux            Colour temperature:4800±200K            Colour reduction index(Ra):96            Total irradiance:310W/m<sup>2</sup>            Size of light field:120~240mm            Rated power of bulb:3.2V/1W            Service life of illuminant:30,000h            240v,50/60Hz         </p>	<p>           Illuminance depth:1050mm            lamp head diameter:500mm            Power supply voltage:AC110-240v,50/60Hz         </p>	<p>           1 piece in a wooden case: 95*92*44 (cm)            GW:37kg            NW:23kg            1 piece in a wooden case: 86*51*24 (cm)            GW:36kg            NW:30kg         </p>
	spring arm germany Shadowless Operating lamp	AC-OL034	<p>           Illuminance:≥150,000Lux/≥150,000Lux            Colour temperature:4500±300K            Colour reduction index(Ra):94            Illuminance depth:≥1000mm            Total irradiance:426W/m<sup>2</sup>            The ratio between irradiance and illuminance :3.3mW/m<sup>2</sup>lx            Size of light field:160~280mm            Rated power of bulb 24V:150W/150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p>	<p>           Colour reduction            Total irradiance:426W/m<sup>2</sup>            Brightness         </p>	<p>           2 piece in a carton: 91*92*37 (cm)            GW:19kg            NW:12kg            1 piece in a wooden case: 122*43*72 (cm)            GW:65kg            NW:50kg         </p>
	spring arm china Shadowless Operating lamp	AC-OL035	<p>           Illuminance:≥150,000Lux /≥150,000Lux            Colour temperature:4500±300K            Colour reduction index(Ra):94            Illuminance depth:≥1000mm            Total irradiance:426W/m<sup>2</sup>            The ratio between irradiance and illuminance :3.3mW/m<sup>2</sup>lx            Size of light field:160~280mm            Rated power of bulb 24V:150W/150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p>	<p>           Colour reduction index(Ra):94            Total irradiance:426W/m<sup>2</sup>            Service life of         </p>	<p>           2 piece in a carton: 91*92*37 (cm)            GW:19kg            NW:12kg            1 piece in a wooden case: 122*43*72 (cm)            GW:65kg            NW:50kg         </p>










	spring arm germany Shadowless Operating lamp	AC-OL036	<p>           Illuminance:≥150,000Lux            Colour temperature:4500±300K            Illuminance depth:≥1000mm            The ratio between irradiance and illuminance :3.3mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p> <p>           Colour reduction index(Ra):94            Total irradiance:426W/m²         </p>	<p>           1 piece in a carton: 91°92*37 (cm)            GW:19kg            NW:12kg            1 piece in a wooden case: 108°43*60 (cm)            GW:46kg            NW:35kg         </p>
	spring arm china Shadowless Operating lamp	AC-OL037	<p>           Illuminance:≥150,000Lux            Colour temperature:4500±300K            Illuminance depth:≥1000mm            The ratio between irradiance and illuminance :3.3mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p> <p>           Colour reduction index(Ra):94            Total irradiance:426W/m²         </p>	<p>           1 piece in a carton: 91°92*37 (cm)            GW:19kg            NW:12kg            1 piece in a wooden case: 108°43*60 (cm)            GW:46kg            NW:35kg         </p>
	TV with camera with monitor (spring arm germany) Shadowless Operating lamp	AC-OL038	<p>           Illuminance:≥150,000Lux            Colour temperature:4500±300K            Colour reduction index(Ra):ZF700:Ra:94/ZF500:Ra:93            Illuminance depth:≥1000mm            Total irradiance: ZF700:426W/m² ZF500:580W/m²            The ratio between irradiance and illuminance :ZF700:3.3mW/m²lx ZF500:4.7mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W/150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic         </p>	<p>           1 piece in a carton: 91°92*37 (cm)            GW:19kg            NW:12kg            1 piece in a carton: 80°81*36 (cm)            GW:14kg NW:8kg            1 piece in a wooden case: 138°85*43 (cm)         </p>
	spring arm germany Shadowless Operating lamp	AC-OL039	<p>           Illuminance:≥150,000Lux            Colour temperature:4500±300K            Colour reduction index(Ra):ZF700:Ra:94/ZF500:Ra:93            Illuminance depth:≥1000mm            Total irradiance: ZF700:426W/m² ZF500:580W/m²            The ratio between irradiance and illuminance :ZF700:3.3mW/m²lx ZF500:4.7mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W/150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic         </p>	<p>           1 piece in a carton: 91°92*37 (cm)            GW:19kg            NW:12kg            1 piece in a carton: 80°81*36 (cm)            GW:14kg NW:8kg            1 piece in a wooden case: 122°43*72 (cm)         </p>
	spring arm china Shadowless Operating lamp	AC-OL040	<p>           Illuminance:≥150,000Lux            Colour temperature:4500±300K            Colour reduction index(Ra):ZF700:Ra:94 ZF500:Ra:93            Illuminance depth:≥1000mm            Total irradiance: ZF700:426W/m² ZF500:580W/m²            The ratio between irradiance and illuminance :ZF700:3.3mW/m²lx ZF500:4.7mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W/150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic         </p>	<p>           1 piece in a carton: 91°92*37 (cm)            GW:19kg            NW:12kg            1 piece in a carton: 80°81*36 (cm)            GW:14kg NW:8kg            1 piece in a wooden case: 122°43*72 (cm)         </p>
	spring arm germany Shadowless Operating lamp	AC-OL041	<p>           Illuminance:≥100,000Lux/≥100,000Lux            Colour temperature:4500±300K            Illuminance depth:≥1000mm            The ratio between irradiance and illuminance :4.7mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W/150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p> <p>           Colour reduction index(Ra):93            Total irradiance:580W/m²         </p>	<p>           2 piece in a carton: 80°81*36 (cm)            GW:14kg NW:8kg            1 piece in a wooden case: 122°43*72 (cm)            GW:65kg            NW:50kg         </p>
	spring arm china Shadowless Operating lamp	AC-OL042	<p>           Illuminance:≥100,000Lux/≥100,000Lux            Colour temperature:4500±300K            Illuminance depth:≥1000mm            The ratio between irradiance and illuminance :4.7mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W/150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p> <p>           Colour reduction index(Ra):93            Total irradiance:580W/m²         </p>	<p>           2 piece in a carton: 80°81*36 (cm)            GW:14kg NW:8kg            1 piece in a wooden case: 122°43*72 (cm)            GW:65kg            NW:50kg         </p>
	spring arm germany Shadowless Operating lamp	AC-OL043	<p>           Illuminance:≥100,000Lux            Colour temperature:4500±300K            Illuminance depth:≥1000mm            The ratio between irradiance and illuminance :4.7mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p> <p>           Colour reduction index(Ra):93            Total irradiance:580W/m²         </p>	<p>           1 piece in a carton: 80°81*36 (cm)            GW:14kg NW:8kg            1 piece in a wooden case: 108°43*60 (cm)            GW:46kg            NW:35kg         </p>
	spring arm china Shadowless Operating lamp	AC-OL044	<p>           Illuminance:≥100,000Lux            Colour temperature:4500±300K            Illuminance depth:≥1000mm            The ratio between irradiance and illuminance :4.7mW/m²lx            Size of light field:160~280mm            Rated power of bulb 24V:150W            Service life of illuminant:1,000h            Power supply voltage:AC110-240v,50/60Hz            Brightness adjustment:Automatic8-stage continuous light adjustment            Shift of major and stand by lamp s:Automatic            Lowest height of installation:2900mm         </p> <p>           Colour reduction index(Ra):93            Total irradiance:580W/m²         </p>	<p>           1 piece in a carton: 80°81*36 (cm)            GW:14kg NW:8kg            1 piece in a wooden case: 108°43*60 (cm)            GW:46kg            NW:35kg         </p>

	Shadowless Operating lamp	AC-OL045	<p>Illuminance:≥120,000Lux/≥120,000Lux  Colour temperature:4500±300K Colour reduction index(Ra):94  Illuminance depth:≥1000mm Total irradiance:426W/m<sup>2</sup>  The ratio between irradiance and illumination :3.3mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W/150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Lowest height of installation:2900mm</p>	<p>2 piece in a carton: 91*92*37 (cm)  GW:19kg  NW:12kg  1 piece in a wooden case: 122*43*72 (cm)  GW:65kg  NW:50kg</p>
	Control Panel Shadowless Operating lamp	AC-OL046	<p>Illuminance:≥120,000Lux  Colour temperature:4500±300K Colour reduction index(Ra):94  Illuminance depth:≥1000mm Total irradiance:426W/m<sup>2</sup>  The ratio between irradiance and illumination :3.3mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Lowest height of installation:2900mm</p>	<p>1 piece in a carton: 91*92*37 (cm)  GW:19kg  NW:12kg  1 piece in a wooden case: 108*43*60 (cm)  GW:46kg  NW:35kg</p>
	Emergency (spring arm germany) Shadowless Operating lamp	AC-OL047	<p>Illuminance:≥100,000Lux  Colour temperature:4500±300K Colour reduction index(Ra):93  Illuminance depth:≥1000mm Total irradiance:580W/m<sup>2</sup>  The ratio between irradiance and illumination :4.7mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Emergency time 2h</p>	<p>1 piece in a wooden case: 95*92*44 (cm)  GW:39kg  NW:26kg  1 piece in a wooden case: 76*53*38 (cm)  GW:44kg  NW:37kg</p>
	Emergency (spring arm china) Shadowless Operating lamp	AC-OL048	<p>Illuminance:≥120,000Lux  Colour temperature:4500±300K Colour reduction index(Ra):93  Illuminance depth:≥1000mm Total irradiance:580W/m<sup>2</sup>  The ratio between irradiance and illumination :4.7mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Emergency time 2h</p>	<p>1 piece in a wooden case: 95*92*44 (cm)  GW:39kg  NW:26kg  1 piece in a wooden case: 76*53*38 (cm)  GW:44kg  NW:37kg</p>
	Emergency Shadowless Operating lamp	AC-OL049	<p>Illuminance:≥120,000Lux  Colour temperature:4500±300K Colour reduction index(Ra):94  Illuminance depth:≥1000mm Total irradiance:426W/m<sup>2</sup>  The ratio between irradiance and illumination :3.3mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Emergency time 2h</p>	<p>1 piece in a wooden case: 95*92*44 (cm)  GW:42kg  NW:29kg  1 piece in a wooden case: 76*53*38 (cm)  GW:44kg  NW:37kg</p>
	Emergency (spring arm germany) Shadowless Operating lamp	AC-OL050	<p>Illuminance:≥150,000Lux  Colour temperature:4500±300K Colour reduction index(Ra):94  Illuminance depth:≥1000mm Total irradiance:426W/m<sup>2</sup>  The ratio between irradiance and illumination :3.3mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Emergency time 2h</p>	<p>1 piece in a wooden case: 95*92*44 (cm)  GW:42kg  NW:29kg  1 piece in a wooden case: 76*53*38 (cm)  GW:44kg  NW:37kg</p>
	Emergency (spring arm china) Shadowless Operating lamp	AC-OL051	<p>Illuminance:≥150,000Lux  Colour temperature:4500±300K Colour reduction index(Ra):94  Illuminance depth:≥1000mm Total irradiance:426W/m<sup>2</sup>  The ratio between irradiance and illumination :3.3mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Emergency time 2h</p>	<p>1 piece in a wooden case: 95*92*44 (cm)  GW:42kg  NW:29kg  1 piece in a wooden case: 76*53*38 (cm)  GW:44kg  NW:37kg</p>
	Emergency Shadowless Operating lamp	AC-OL052	<p>Illuminance:≥100,000Lux Colour temperature:4500±300K  Colour reduction index(Ra):93 Illuminance depth:≥1000mm  Total irradiance:580W/m<sup>2</sup>  The ratio between irradiance and illumination :4.7mW/m<sup>2</sup>lx  Size of light field:160~280mm  Rated power of bulb 24V:150W Service life of illuminant:1,000h  Power supply voltage:AC110-240v,50/60Hz  Brightness adjustment:Automatic8-stage continuous light adjustment  Shift of major and stand by lamp s:Automatic  Lowest height of installation:2900mm  Emergency time 2h</p>	<p>1 piece in a carton: 91*92*37 (cm)  GW:19kg  NW:12kg  1 piece in a wooden case: 108*43*60 (cm)  GW:46kg  NW:35kg  1 piece in a wooden case:</p>
	Cold light operating lamp	AC-OL053	<p>Illuminance:5-reflector:≥50,000Lux  12-reflector:≥120,000Lux  Colour temperature:4000±500K  Size of light field:160~220mm  Illuminance depth:≥600mm  Total irradiance:5-reflector:490W/m<sup>2</sup> 12-reflector:718W/m<sup>2</sup>  Rated power of bulb:24V/25W  Power supply voltage:220V/50Hz  Rated power:450W  Lowest height of installation:2900mm  Bulb:Cold light source</p>	<p>1 piece in a wooden case:70*70*38 (cm)  GW:17kg NW:8kg  1 piece in a wooden case:114*114*43 (cm)  GW:40kg  NW:22kg  1 piece in a</p>








	Cold light operating lamp	AC-OL054	<p>Illuminance:6-reflector:60,000Lux~120,000Lux 3-reflector:40,000Lux~80,000Lux Colour temperature:4000±500K Size of light field:160~220mm Illuminance depth:≥600mm Total irradiance:6-reflector:718W/m<sup>2</sup> 3-reflector:490W/m<sup>2</sup> Rated power of bulb:24V/50W Power supply voltage:220V/50Hz Rated power:660W Lowest height of installation:2900mm Bulb:Tungsten halogen bulb</p>	<p>1 piece in a wooden case:70*70*38 (cm) GW:17kg NW:8kg 1 piece in a wooden case:114*114*43 (cm) GW:40kg NW:22kg 1 piece in a</p>
	Cold light operating lamp	AC-OL055	<p>Illuminance:≥50,000Lux/≥50,000Lux Colour temperature:4000±500K Size of light field:160~220mm Illuminance depth:≥600mm Total irradiance:490W/m<sup>2</sup> Rated power of bulb:24V/25W Power supply voltage:220V/50Hz Rated power:275W Lowest height of installation:2900mm Bulb:Cold light source</p>	<p>1 piece in a wooden case:70*70*38 (cm) GW:17kg NW:8kg 1 piece in a wooden case:114*43*72 (cm) GW:78kg NW:6kg</p>
	Cold light operating lamp	AC-OL056	<p>Illuminance:≥120,000Lux Colour temperature:4000±500K Size of light field:160~220mm Illuminance depth:≥600mm Total irradiance:718W/m<sup>2</sup> Rated power of bulb:24V/25W Power supply voltage:220V/50Hz Rated power:350W Lowest height of installation:2900mm Bulb:Cold light source</p>	<p>1 piece in a wooden case:114*114*43 (cm) GW:40kg NW:22kg 1 piece in a wooden case:104*43*57 (cm) GW:48kg NW:38kg</p>
	Cold light operating lamp	AC-OL057	<p>Illuminance:≥90,000Lux Colour temperature:4000±500K Size of light field:160~220mm Illuminance depth:≥600mm Total irradiance:718W/m<sup>2</sup> Rated power of bulb:24V/25W Power supply voltage:220V/50Hz Rated power:240W Lowest height of installation:2900mm Bulb:Cold light source</p>	<p>1 piece in a wooden case:114*114*43 (cm) GW:38kg NW:21kg 1 piece in a wooden case:104*43*57 (cm) GW:48kg NW:38kg</p>
	Cold light operating lamp	AC-OL058	<p>Illuminance:60,000Lux~120,000Lux Colour temperature:4000±500K Size of light field:160~220mm Illuminance depth:≥600mm Total irradiance:718W/m<sup>2</sup> Rated power of bulb:24V/50W Power supply voltage:220V/50Hz Rated power:350W Lowest height of installation:2900mm Bulb:Tungsten halogen bulb</p>	<p>1 piece in a wooden case:114*114*43 (cm) GW:38kg NW:21kg 1 piece in a wooden case:104*43*57 (cm) GW:48kg NW:38kg</p>
	Cold light operating lamp	AC-OL059	<p>Illuminance:≥50,000Lux Colour temperature:4000±500K Size of light field:160~220mm Illuminance depth:≥600mm Total irradiance:490W/m<sup>2</sup> Rated power of bulb:24V/25W Power supply voltage:220V/50Hz Rated power:150W Lowest height of installation:2900mm Bulb:Cold light source</p>	<p>Meas:81*69*33 (cm) GW:30kg NW:19kg 1 piece in a wooden case</p>
	Cold light operating lamp	AC-OL060	<p>Illuminance:≥50,000Lux Colour temperature:4000±500K Colour temperature:160~220mm Illuminance depth:≥600mm Total irradiance:490W/m<sup>2</sup> Rated power of bulb:24V/25W Power supply voltage:220V/50Hz Rated power:150W Lowest height of installation:2300-2600mm Bulb:Cold light source</p>	<p>Meas:90*78*34 (cm) GW:41kg NW:29kg 1 piece in a wooden case</p>
	Cold light operating lamp	AC-OL061	<p>Illuminance:≥40,000Lux Colour temperature:4000±500K Colour temperature:160~220mm Illuminance depth:≥600mm Total irradiance:490W/m<sup>2</sup> Rated power of bulb:24V/25W Power supply voltage:220V/50Hz Rated power:120W Lowest height of installation:2900mm Bulb:Cold light source</p>	<p>Meas:81*69*33 (cm) GW:30kg NW:19kg 1 piece in a wooden case</p>
	Cold light operating lamp	AC-OL062	<p>Illuminance:40,000Lux~80,000Lux Colour temperature:4000±500K Colour temperature:160~220mm Illuminance depth:≥600mm Total irradiance:490W/m<sup>2</sup> Rated power of bulb:24V/50W Power supply voltage:220V/50Hz Rated power:200W Lowest height of installation:2900mm Bulb:Tungsten halogen bulb</p>	<p>Meas:81*69*33 (cm) GW:30kg NW:19kg 1 piece in a wooden case</p>

	Control Panel Cold light operating lamp	AC-OL063	<p>Illuminance:40,000Lux~80,000Lux  Colour temperature:4000±500K  Colour temperature:160 ~ 220mm  Illuminance depth:≥600mm  Total irradianc:490W/m<sup>2</sup>  Rated power of bulb:24V/50W  Power supply voltage:220V/50Hz  Rated power:200W  Lowest height of installation:2900mm  Bulb:Tungsten halogen bulb</p>	<p>Meas:100*85*38  (cm)  GW:45kg  NW:38kg  1 piece in a wooden case</p>
	on wall (Control Panel) Cold light operating lamp	AC-OL064	<p>Illuminance:40,000Lux~80,000Lux  Colour temperature:4000±500K  Colour temperature:160 ~ 220mm  Illuminance depth:≥600mm  Total irradianc:490W/m<sup>2</sup>  Rated power of bulb:24V/50W  Power supply voltage:220V/50Hz  Rated power:200W  Lowest height of installation:2300-2600mm  Bulb:Tungsten halogen bulb</p>	<p>Meas:110*84*43  (cm)  GW:50kg  NW:37kg  1 piece in a wooden case</p>
	Cold light operating lamp	AC-OL065	<p>Illuminance :≥60,000Lux  color temperature :3700±500K  Illuminance depth :≥600mm  Total irradiance :547W/m<sup>2</sup>  size of light field:160-220mm  Rated power of bulb:25W/24v  import power :150VA  Power supply voltage:220V/50Hz  bulb :cold light</p>	<p>Meas:114*64*33  (cm)  GW:54kgNW:42kg  1 piece in a wooden case</p>
	Cold light operating lamp	AC-OL066	<p>Illuminance :≥50,000Lux  color temperature :3700±500K  Illuminance depth :≥600mm  Total irradiance :547W/m<sup>2</sup>  size of light field:160-220mm  Rated power of bulb:25W/24v  Import power :150VA  Power supply voltage:220V/50Hz  bulb :cold light</p>	<p>Meas:114*64*33  (cm)  GW:54kg  MW:42kg  1 piece in a wooden case</p>
	Cold light operating lamp	AC-OL067	<p>Illuminance:40,000~80,000Lux  color temperature:3700±500K  Illuminance dept:≥600mm  Total irradiance:547W/m<sup>2</sup>  size of light field:160-220mm  Rated power of bulb:24V/50W  Import power:200VA  Power supply voltage:220V/50Hz  bulb :cold light</p>	<p>Meas:114*64*33  (cm)  GW:54kg  NW:42kg  1 piece in a wooden case</p>
	Emergency Cold light operating lamp	AC-OL068	<p>Illuminance :≥60,000Lux  color temperature :3700±500K  Illuminance depth :≥600mm  Total irradiance :547W/m<sup>2</sup>  size of light field:160-220mm  Rated power of bulb:25W/24v  import power :150VA  Power supply voltage:220V/50Hz  bulb :cold light</p> <p>Emergency time:3 hours</p>	<p>Meas:114*65*38  (cm)  GW:65kg  NW:52kg  1 piece in a wooden case</p>
	Emergency Cold light operating lamp	AC-OL069	<p>Illuminance :≥50,000Lux  color temperature :3700±500K  Illuminance depth :≥600mm  Total irradiance :547W/m<sup>2</sup>  size of light field:160-220mm  Rated power of bulb:25W/24v  import power :150VA  Power supply voltage:220V/50Hz  bulb :cold light</p> <p>Emergency time:3 hours</p>	<p>Meas:114*65*38  (cm)  GW:65kg  NW:52kg  1 piece in a wooden case</p>
	Emergency Cold light operating lamp	AC-OL070	<p>Illuminance :≥50,000Lux  color temperature :3700±500K  Illuminance depth :≥600mm  Total irradiance :547W/m<sup>2</sup>  size of light field:160-220mm  Rated power of bulb:25W/24v  Import power :150VA  Power supply voltage:220V/50Hz  bulb :cold light</p> <p>Emergency time:3 hours</p>	<p>Meas:114*65*38  (cm)  GW:65kg  NW:52kg  1 piece in a wooden case</p>
	Adjust light Cold light operating lamp	AC-OL071	<p>Illuminance :≥25000Lux  color temperature :3500±500K  Power supply voltage:220V/50Hz  Rated power of bulb:24V/50W  Import power:56VA  bulb :Halogen bulb</p>	<p>Meas:124*27*27  (cm)  GW:12kg NW:9kg  1 piece in a carton</p>



	Cold light operating lamp	AC-OL072	<p>Illuminance :<math>\geq</math>12000Lux                      color temperature :3500<math>\pm</math>500K  Power supply voltage:220V/50Hz  Rated power of bulb:24V/25W  Import power:56VA  bulb :cold light bulb</p>	<p>Meas:124*27*27  (cm)  GW:12kg NW:9kg  1 piece in a wooden case</p>
	(Emergency) Cold light operating lamp	AC-OL073	<p>Illuminance :<math>\geq</math>12000Lux                      color temperature :3500<math>\pm</math>500K  Power supply voltage:220V/50Hz  Rated power of bulb:24V/25W  Import power:56VA  bulb :cold light bulb    Emergency time:3 hours</p>	<p>Meas:123*33*23  (cm)  GW:19kg  NW:14kg  1 piece in a wooden case</p>
	(On ceiling) cold light Operating lamp (deep)	AC-OL074	<p>Illuminance:<math>\geq</math>50,000Lux  color temperature :4000<math>\pm</math>500k  Power supply voltage:220V/50Hz  Rated power of bulb :24V/50W  Bulb:halogen bulb  Lowest height of installation :2300-2600mm</p>	<p>Meas:110*34*41  (cm)  GW:41kg  NW:29kg  1 piece in a wooden case</p>
	(On wall) cold light Operating lamp (deep)	AC-OL075	<p>Illuminance:<math>\geq</math>50,000Lux  color temperature :4000<math>\pm</math>500k  Power supply voltage :220V/50Hz  Rated power of bulb :24V/50W  Bulb:halogen bulb  Lowest height of installation :2300-2600mm</p>	<p>Meas:120*52*30  (cm)  GW:36kg  NW:28kg  1 piece in a wooden case</p>
	Cold light Operating lamp (deep)	AC-OL076	<p>Illuminance:<math>\geq</math>50,000Lux  color temperature :4000<math>\pm</math>500k  Power supply voltage:220V/50Hz  Rated power of bulb :24V/50W  Bulb:halogen bulb</p>	<p>Meas:118*51*29  (cm)  GW:45kg  NW:40kg  1 piece in a wooden case</p>
	Emergency cold light Operating lamp (deep)	AC-OL077	<p>Illuminance:<math>\geq</math>50,000Lux  color temperature :4000<math>\pm</math>500k  Power supply voltage:220V/50Hz  Rated power of bulb :24V/50W  Bulb:halogen bulb                      Emergency time:3 hours</p>	<p>Meas:120*54*38  (cm)  GW:56kg  NW:43kg  1 piece in a wooden case</p>
	Examination lamp	AC-OL078	<p>Illuminance :20W:20,000Lux 0.5m/5,000Lux 1m  35W:25,000Lux 0.5m/10,000Lux 1m  50W:30,000Lux 0.5m/20,000Lux 1m  Color temperature :4000<math>\pm</math>200K  Size of light field :100mm  Adjust : adjust light  Power supply voltage :220V/50Hz  Rated power of bulb :12V/20W,12V/35W,12V/50W  crutch:adjust up/down  bulb: halogen bulb</p>	<p>Meas:95*56*17  (cm)  GW:6.5kg  NW:5.2kg  1 piece in a wooden case</p>
	Examination lamp	AC-OL079	<p>Illuminance :20W:20,000Lux 0.5m/5,000Lux 1m  35W:25,000Lux 0.5m/10,000Lux 1m  50W:30,000Lux 0.5m/20,000Lux 1m  Color temperature :4000<math>\pm</math>200K  Size of light field :100mm  Adjust : adjust light  Power supply voltage :220V/50Hz  Rated power of bulb :12V/20W,12V/35W,12V/50W  bulb: halogen bulb</p>	<p>Meas:75*19*25  (cm)  GW:2.9kg  NW:2.1kg  1 piece in a carton  Meas :75*19*25  (cm) GW:5kg  NW:4.2kg  2 piece in a carton</p>
	LED cold light Operating lamp (deep)	AC-OL080	<p>Illuminance:<math>\geq</math>25,000Lux  color temperature :4800<math>\pm</math>500k  Power supply voltage:220V/50Hz  Rated power of bulb :3.2V/1W  Bulb:LED bulb</p>	<p>Meas:118*51*29  (cm)  GW:45kg  MW:40kg  1 piece in a carton</p>

	LED (Emergency)cold light Operating lamp (deep)	AC-OL081	<p>Illuminance:≥25,000Lux color temperature :4800±500K Power supply voltage:220V/50Hz Rated power of bulb :3.2V/1W Bulb:LED bulb</p> <p>Emergency time:3 hours</p>	<p>Meas:120*54*38 ( cm) GW:56kg NW:43kg 1 piece in a wooden case</p>
	LED Examination lamp	AC-OL082	<p>Illuminance :30,000Lux 0.5m/5,000Lux 1m Color temperature :4800±200K Size of light field :100mm Adjust : adjust light Power supply voltage :220V/50Hz Rated power of bulb :3.3/3w Crutch:adjust up/down Bulb:LED bulb</p>	<p>Meas:95*56*17 ( cm) GW:6.5kg NW:5.2kg 1 piece in a carton</p>
	LED (on wall) Examination lamp	AC-OL083	<p>Illuminance :30,000Lux 0.5m/5,000Lux 1m Color temperature :4800±200K Size of light field :100mm Adjust : adjust light Power supply voltage :220V/50Hz Rated power of bulb :3.3/3w Bulb:LED bulb</p>	<p>Meas:75*19*25 ( cm) GW:2.9kg NW:2.1kg 1 piece in a carton Meas:75*19*25 ( cm) GW:5kg NW:4.2kg 2 piece in a carton</p>
	Operating lamp LED	AC-OL084	<p>Illuminance:≥ 20,000Lux color temperature:4800±200K Size of light field:120mm Power supply voltage:220V/50Hz Rated power of bulb:3.2V1W lamp head diameter:160mm</p>	<p>Meas:75*20*33 ( cm) GW:5.3kg NW:4.6kg 1 piece in a carton</p>
	LED Operating lamp	AC-OL085	<p>Illuminance:≥ 20,000Lux color temperature:4800±200K Size of light field:120mm Power supply voltage:220V/50Hz Rated power of bulb:3.2V1W lamp head diameter:160mm Emergency time:5 hours</p>	<p>Meas:75*20*33 ( cm) GW:5.3kg NW:4.6kg 1 piece in a carton</p>
	LED (on ceiling) cold light Operating lamp (deep)	AC-OL086	<p>Illuminance:≥25,000Lux color temperature :4800±500K Power supply voltage:220V/50Hz Rated power of bulb :3.2V/1wW Bulb:LED bulb Lowest height of installation :2300-2600mm</p>	<p>Meas:110*34*41 ( cm) GW:41kg NW:29kg 1 piece in a wooden case</p>
	LED (On wall) cold light Operating lamp (deep)	AC-OL087	<p>Illuminance:≥25,000Lux color temperature :4800±500K Power supply voltage:220V/50Hz Rated power of bulb :3.2V/1wW Bulb:LED bulb Lowest height of installation :2300-2600mm</p>	<p>Meas:120*52*30 ( cm) GW:36kg NG:28kg 1 piece in a wooden case</p>
	LED cold light Operating lamp	AC-OL088	<p>Illuminance :≥25000Lux color temperature :4800±500K Power supply voltage:220V/50Hz Rated power of bulb:3.2V/1W Import power:56VA bulb :LED bulb LED bulb life:≥=30000h color reduction index(Ra): 93</p>	<p>Meas:124*27*27 ( cm) GW:12kg NW:9kg 1 piece in a carton</p>
	LED Emergency cold light Operating lamp	AC-OL089	<p>Illuminance :≥25000Lux color temperature :4800±500K Power supply voltage:220V/50Hz Rated power of bulb:3.2V/1W Import power:56VA bulb :LED bulb LED bulb life:≥=30000h color reduction index(Ra): 93 Emergency time:3 hours</p>	<p>Meas:123*33*23 ( cm) GW:19kg NW:14kg 1 piece in a wooden case</p>

	LED cold light Operating lamp	AC-OL090	<p>Illuminance(Lux):≥80,000  color temperature(K):4800±300  illumination depth (mm):≥500  Size of light field(mm):160-220  Rated power of bulb :3.2V/1W (32pc led bulb )  Power supply voltage:~220V/50Hz  LED bulb life (hr) :≥30,000</p>	<p>Meas:114*64*33  (cm)  GW:54kg  NW:42kg  1 piece in a  wooden case</p>
	LED (Common Arm) cold light Operating lamp	AC-OL091	<p>Illuminance(Lux):≥100,000  color temperature(K):4800±300  illumination depth (mm):≥500  Size of light field(mm):160-220  Rated power of bulb :3.2V/1W (32pc led bulb )  Power supply voltage:~220V/50Hz  LED bulb life (hr) :≥30,000</p>	<p>Meas:90*69*32  (cm)  GW:45kg  NW:34kg  1 piece in a  wooden case</p>
	(on wall) cold light Operating lamp	AC-OL092	<p>Illuminance:40,000~80,000Lux  color temperature:3700±500K  illumination dept:≥600mm  Total irradiance:547W/m²  size of light field:160-220mm  Rated power of bulb:24V/50W  Import power:200VA  Power supply voltage:220V/50Hz  bulb :cold light Lowest height of installation :2300-2600mm</p>	<p>Meas:110*84*43  (cm)  GW:50kg  NW:37kg  1 piece in a  wooden case</p>
	LED (Common Arm) cold light Operating lamp	AC-OL093	<p>Illuminance(Lux):≥80,000  color temperature(K):4800±300  illumination depth (mm):≥500  Size of light field(mm):160-220  Rated power of bulb :3.2V/1W (32pc led bulb )  Power supply voltage:~220V/50Hz  LED bulb life (hr) :≥30,000</p>	<p>Meas:114*64*33  (cm)  GW:54kg  NW:42kg  1 piece in a  wooden case</p>
	LED (Common Arm) cold light Operating lamp	AC-OL094	<p>Illuminance(Lux):≥100,000  color temperature(K):4800±300  illumination depth (mm):≥500  Size of light field(mm):160-220  Rated power of bulb :3.2V/1W (32pc led bulb )  Power supply voltage:~220V/50Hz  LED bulb life (hr) :≥30,000</p>	<p>Meas:114*64*33  (cm)  GW:54kg  NW:42kg  1 piece in a  wooden case</p>
	LED (Common Arm) cold light Operating lamp	AC-OL095	<p>Illuminance(Lux):≥100,000  color temperature(K):4800±300  illumination depth (mm):≥500  Size of light field(mm):160-220  Rated power of bulb :3.2V/1W (32pc led bulb )  Power supply voltage:~220V/50Hz  LED bulb life (hr) :≥30,000 Total irradiance 480W/M2  height of installation :2300-2600mm</p>	<p>Meas:90*69*32  (cm)  GW:45kg  NW:34kg  1 piece in a  wooden case</p>
	LED (Common Arm) cold light Operating lamp	AC-OL096	<p>Illuminance(Lux):≥100,000/100,000  color temperature(K):4800±300  illumination depth (mm):≥500  Size of light field(mm):160-220  Rated power of bulb :3.2V/1W (32pc led bulb )  Power supply voltage:~220V/50Hz  LED bulb life (hr) :≥30,000</p>	<p>lamp head  :63*63*32 (cm)  GW:14kg NW:8kg  2 piece in a  wooden case  RAM:90*35*58  (cm)  GW:46kg  NW:36kg  1 piece in a  wooden case</p>