










ZHANGJIAGANG ANNECY MACHINERY CO., LTD.

**Stainless Steel Trolley Catalogue & Quotation**

Images	Model No	Details
	AC-SPT001	<p>Specification: 730*480*800mm</p> <p>1. The tabletop is formed by ABS injection molding process at one time; double-layer tabletop, three-layer stainless steel double-layer guardrail for each layer of tabletop, four-pillar stainless steel material 1.0 thickness, Φ25mm, wear-resistant, easy to clean, 120mm in one middle, inner space: 430*335* 110mm, 3*3 separators in the drawer, can be separated freely, the drawer handle is dovetail type, sealed slot type transparent identification card Specifications: 115*28mm, to prevent liquid and dust from entering;</p> <p>2. The bottom of the car body: universal mute wheels, two of which have brake function.</p> <p>3. Packing size: 860*310*500mm</p> <p>4. Gross weight: 12.75kg</p>
	AC-SPT002	<p>Specification: 730*480*985mm</p> <p>1. The table top is formed by ABS injection process at one time; three-layer table top, three-layer stainless steel double-layer guardrail for each layer, four-post stainless steel material 1.0 thickness, Φ25mm, wear-resistant and easy to clean;</p> <p>2. The bottom of the car body: universal mute wheels, two of which have brake function.</p> <p>3. Packing size: 860*310*500mm</p> <p>4. Gross weight: 14.3kg</p>
	AC-SPT003	<p>Specification: 730*480*985mm</p> <p>1. The tabletop is formed by ABS injection molding process at one time; double-layer tabletop, three-layer stainless steel double-layer guardrail for each layer of tabletop, four-pillar stainless steel material 1.0 thickness, Φ25mm, wear-resistant, easy to clean, Specifications: 115*28mm, to prevent liquid and dust from entering;</p> <p>2. The bottom of the car body: universal mute wheels, two of which have brake function.</p> <p>3. Packing size: 860*310*500mm</p> <p>4. Gross weight: 12.75kg</p>
	AC-SPT004	<p>1. The tabletop is formed by ABS injection molding process at one time; double-layer tabletop, three-layer stainless steel double-layer guardrail for each layer of tabletop, four-pillar stainless steel material 1.0 thickness, Φ25mm, wear-resistant, easy to clean, 120mm in one middle, inner space: 430*335* 110mm, 3*3 separators in the drawer, can be separated freely, the drawer handle is dovetail type, sealed slot type transparent identification card Specifications: 115*28mm, to prevent liquid and dust from entering;</p> <p>2. The bottom of the car body: universal mute wheels, two of which have brake function. And with one bucket</p> <p>3. Packing size: 860*310*500mm</p> <p>4. Gross weight: 12.75kg</p>
	AC-SPT005	<p>1. The tabletop is formed by ABS injection molding process at one time; double-layer tabletop, three-layer stainless steel double-layer guardrail for each layer of tabletop, four-pillar stainless steel material 1.0 thickness, Φ25mm, wear-resistant, easy to clean, 120mm in one middle, inner space: 430*335* 110mm, 3*3 separators in the drawer, can be separated freely, the drawer handle is dovetail type, sealed slot type transparent identification card Specifications: 115*28mm, to prevent liquid and dust from entering;</p> <p>2. The bottom of the car body: universal mute wheels, two of which have brake function. And with one bucket and one stainless steel pot</p> <p>3. Packing size: 860*310*500mm</p> <p>4. Gross weight: 12.75kg</p>
	AC-SPT006	<p>Specification: 640*465*995mm</p> <p>1. The tabletop is formed by ABS injection molding process at one time; double-layer tabletop, three-layer stainless steel double-layer guardrail for each layer of tabletop, four-pillar stainless steel material 1.0 thickness, Φ25mm, wear-resistant, easy to clean, 120mm in one middle, inner space: 430*335* 110mm, 3*3 separators in the drawer, can be separated freely, the drawer handle is dovetail type, sealed slot type transparent identification card Specifications: 115*28mm, to prevent liquid and dust from entering;</p> <p>2. The bottom of the car body: universal mute wheels, two of which have brake function. And with two buckets</p> <p>3. Packing size: 860*310*500mm</p> <p>4. Gross weight: 12.75kg</p>
	AC-SPT007	<p>Specification: 640*465*995mm</p> <p>1. The tabletop is formed by ABS injection molding process at one time; double-layer tabletop, three-layer stainless steel double-layer guardrail for each layer of tabletop, four-pillar stainless steel material 1.0 thickness, Φ25mm, wear-resistant, easy to clean, 120mm in one middle, inner space: 430*335* 110mm, 3*3 separators in the drawer, can be separated freely, the drawer handle is dovetail type, sealed slot type transparent identification card Specifications: 115*28mm, to prevent liquid and dust from entering;</p> <p>2. The bottom of the car body: universal mute wheels, two of which have brake function.</p> <p>3. Packing size: 860*310*500mm</p> <p>4. Gross weight: 12.75kg</p>