













ZHANGJIAGANG ANNECY MACHINERY CO., LTD.

CATALOGUE & QUOTATION

Images	Description	Model No	Details
	Medical pendulum	AC-PD001	Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 10 pcs Instrument platform: 2 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request
	Medical pendulum	AC-PD002	Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 10 pcs Instrument platform: 2 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request
	Medical pendulum	AC-PD003	Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 10 pcs Instrument platform: 3 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request
	Medical pendulum	AC-PD004	Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 10 pcs Instrument platform: 3 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request
	Medical pendulum	AC-PD005	Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 6 pcs Instrument platform: 2 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request

	<p>Medical pendulum</p>	<p>AC-PD006</p>	<p>Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 6 pcs Instrument platform: 2 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request</p>
	<p>Medical pendulum</p>	<p>AC-PD007</p>	<p>Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 10 pcs Instrument platform: 3 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request</p>
	<p>Medical pendulum</p>	<p>AC-PD008</p>	<p>Gas outlet Dry side: Oxygen*2, Vacuum*2, Compressed Air*2; Wet side: Oxygen*2, Vacuum*2, Compressed Air*2; Infusion Support Dry side: 1 set; Wet side: 1 set; Instrument platform Dry side: 2 layers; Wet side: 3 layers; Grounding terminal Dry side: 2 layers; Wet side: 2 layers; Electrical sockets: 10 pcs; Network and communication interface: 1 each Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request</p>
	<p>Medical pendulum</p>	<p>AC-PD009</p>	<p>Gas Outlet for Oxygen: 2* Gas Outlet for Vacuum: 2* Gas Outlet for Compressed Air: 2* Infusion Support: 1 set Electrical sockets: 10 pcs Instrument platform: 3 layers Network and communication interface: 1 each Grounding terminal: 2 pcs Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request</p>
	<p>Medical pendulum</p>	<p>AC-PD010</p>	<p>Gas outlet Dry side: Oxygen*2, Vacuum*2, Compressed Air*2; Wet side: Oxygen*2, Vacuum*2, Compressed Air*2; Infusion Support Dry side: 1 set; Wet side: 1 set; Instrument platform Dry side: 2 layers; Wet side: 3 layers; Grounding terminal Dry side: 2 layers; Wet side: 2 layers; Electrical sockets: 10 pcs; Network and communication interface: 1 each Medical Gas Outlets Standards: Optional (British, German, French, Metric, Ohmeda, DISS, Chemetron, Purit-Bennett) Electrical Sockets: Customized as clients' request</p>